Hinri-Alloy CB

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Hinri-Alloy CB

Instructions for use

Area of application: Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 4, for dental devices.

Product description: Hinri-Alloy CB is a non-precious, corrosion resistant, cobalt-based bonding alloy. Hinri-Alloy CB is applicable for the laser technique and free of nickel and beryllium.

Composition (in % by mass):

Co: 61.5 Cr: 28.0 W 8.5 Si: 1.7

Other constituents: Fe, Mn, C

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines):

Proof stress 0.2 % (MPa):	620
Elongation at rupture (%):	10
Modulus of elasticity (GPa):	190
Coefficient of expansion:	20°C-600°C $14.1 \times 10^{-6}/K$
Tensile strength (MPa):	845
Vickers hardness VH 10:	285
Density (g/cm ³):	8.3
Melting range (°C):	1.390-1.415
Max. oxide firing temperature (°C):	980
Casting temperature (°C):	1.515

Investing and casting: Hinri-Alloy CB can be used with phosphate-bonded crown and bridge investments, such as Hinrivest[®] KB. It is preheated to 850 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the ERNST HINRICHS Technique instructions when attaching sprues. Use a separate ceramic crucible for casting Hinri-Alloy CB high-grade, non-precious alloy to prevent it being contaminated by other alloys. Clean the crucible after each cast. Cylindrical, continuous cast ingots collapse more easily during induction melting of the alloy. The alloy is cast when all the cylindrical ingots have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. Please note that the melting point may vary depending on the type of casting machine used and when melting with vacuum, for example, less oxide forms and the oxide layer in the vacuum disintegrates more quickly. Considerable experience is required when melting high-grade alloys automatically, and it is most important that heating is set exactly according to the manufacturer's instructions. Inaccurate heat setting impairs the alloy and can have a detrimental effect on the porcelain bond. Always use the reducing zone of the flame in a circular motion to melt the alloy. Casting Hinri-Alloy CB repeatedly is not recommended, as this could alter both its composition and properties.

Firing porcelain: The minimum thickness of alloy before firing porcelain is 0.2 mm. Before applying the ceramic, fire at 960 – 980 °C, sandblast with 110 µm aluminum oxide and clean carefully with steamer. **Attention:** Never put Hinri-Alloy CB in acid bath. Long term cooling after firing is generally recommended.

Soldering / welding: For soldering of Hinri-Alloy CB all commercial cobalt based solders are recommended.

Finishing / Polishing: Thru its low hardness of Hinri-Alloy CB the surface is easy to finish and polish.

Safety hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask (Recommendation: type FFP2)!

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Storage conditions: Dry, at normal conditions.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

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Commerciale: ERNST HINRICHHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH segnare il confezionamento a **CE 0124**
ERNST HINRICHHS GmbH