Hinri-Alloy CoCr

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Hinri-Alloy CoCr

Instructions for use

Area of application: Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5

Product description: Hinri-Alloy CoCr is an alloy for denture bases is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, Moreover its molecular structure allows to obtain smooth, compact surfaces with little oxide formation. It is highly corrosion resistant and is free of toxic elements like beryllium, indium and gallium.

Composition (in % by mass):

Co: 62.0 Cr: 31.0 Mo 5.0

Other constituents: Si, C, Mn

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines):

Proof stress 0,2 % (MPa):	650
Elongation at rupture (%)	5.0
Modulus of elasticity (GPa):	220
Tensile strength (MPa):	890
Vickers hardness VH 10:	350
Density (g/cm ³):	8.3
Melting range (°C):	1.280 - 1.360
Casting temperature (°C):	1.340

Investing and casting: Hinri-Alloy CoCr can be used with phosphate-bonded CrCo investments such as Hinriinvest[®] CoCr or Hinriinvest[®] SG. Prepare the moulds for casting according to the ERNST HINRICHS ring less system. Preheating range is 850 °C – 1.000 °C. We recommend 920 °C. Heat soak at casting temperature for 30 minutes. Use a separate ceramic crucible for casting Hinri-Alloy CoCr high-grade, non-precious alloy to prevent it being contaminated by other alloys. Clean the crucible after each cast.

The alloy is cast with an induction machine when all the cylindrical ingots have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. Please note that the melting point may vary depending on the type of casting machine used and when melting with a vacuum, for example, less oxide forms and the oxide layer in the vacuum disintegrates more quickly. Considerable experience is required when melting high-grade alloys automatically, and it is most important that the heating is set exactly according to the manufacturer's instructions. Accurate heat setting is essential to protect the properties of the alloy. Light the flame and adjust the pressure in order to have a blue nucleus about 4 mm long. Place some metal ingots into the preheated crucible. Move the flame over the crucible rotating the tip, in order to evenly heat the alloy. As soon as the ingots are wholly melted, and after disappearing of their shadows, start centrifugal action. Casting Hinri-Alloy CoCr repeatedly is not recommended. During re-use pay attention that the alloys are from the same batch.

Preparing and polishing: Cool down the muffle to room temperature, Sandblast the framework with aluminum oxide 120-150 µm and continue the working cycle according to your personal method. Remove carefully any small casting blows, go over all the corners and edges with a round bur and then smooth the surfaces with rubber abrasives and polish them.

Soldering and welding: Presolder using CrCo solder (2 mm, REF 893060, 1 mm, REF 893061) and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys. Laser weld using a CrCo laser (REF 893062) welding rod.

Safety hint: Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask.

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

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Commerciale: ERNST HINRICHHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH segnare il confezionamento a **CE 0124**
ERNST HINRICHHS GmbH