1. Usage

US Splint Memory milling blanks are used to manufacture occlusal splints with thermomemory effect.

2. Product description

The Splint Memory material is characterized by its thermoplastic flexibility resulting in highly precise adaption to the tooth situation. Moreover, Splint Memory offers stress-free wearing comfort for the patient. Furthermore, the self-adjusting occlusal splint shows high transparency and is also extremely resilient to breakage. The industrial production process of Splint Memory blanks provides an exceptional material homogeneity and thereby an outstanding long-term stability. The usage of Splint Memory milling blanks by CAD/CAM technology ensures a safe production process. The use of Splint Memory milling blanks using CAD/CAM technology also ensures a safe process, since mixing errors (e.g. by hand mixing) are eliminated. In addition the odor occurring during mixing will be avoided. Indications: Material for the production of dental occlusal splints: Reflex splints, Therapeutic splints, Repositioning splints, Stabilization splints.

3. Processing

■ Splint Memory milling blanks should be used only for the purpose of CAD/CAM production of dental occlusal splints. They can be used in all common CAD/CAM systems (**fig. 1**).
 ■ Suitable cross-cut carbide cutters are to be used exclusively for the processing of the Splint Memory material.
 ■ Splint Memory occlusal splints are designed and milled by CAD/CAM technology and qualified staff. After finishing, the thickness in the occlusal region of the splint must not be less than 0,9 mm. For aesthetic reasons, a labial reduction to 0,8 mm is permissible. Any corners and edges should be rounded off.
 ■ Splint Memory splints can be removed from the milling blank by using suitable cross-cut carbide cutters or cutting discs suitable for composites. Carefully cut through the retaining strips without pressure.
 ■ In order to avoid plaque accumulation subsequent polishing of the splints is essential. Please pre-polish with suitable silicone polishers and goat's hair brushes. The high-luster polishing should be done with a corresponding polishing compound (**fig. 2**). As far as possible, please avoid heat generation during the polishing and finishing of the Splint Memory splint. This guarantees the optimal fitting of the splint.
 ■ **Please note:** Upon completion, use an ultrasonic bath with water at a maximum temperature of approx. 35 °C for cleaning. Cleaning solutions are not recommended.
 ■ Splint Memory powder liquid system must be used for repairs or addition of material in order to preserve the thermoplastic properties of the splint. The processing recommendations of the manufacturer must be taken in consideration.

4. Milling parameters Splint Memory

Roughing:

- 2 mm ball radius cutter – single cutter
- Speed: 22,000 rpm
- Radial path distance: 0.6 mm
- Z depth increase: 0.6 mm
- Feed: 1,600 mm/min

Finishing:

- 2 mm ball radius cutter – single cutter
- Speed: 22,000 rpm
- Radial path distance: 0.1 mm
- Z depth increase: 0.1 mm
- Feed: 1,600 mm/min

Rest roughing:

- 1 mm ball radius cutter – single cutter
- Speed: 28,000 rpm
- Radial path distance: 0.1 mm
- Z depth increase: 0.1 mm
- Feed: 1,200 mm/min


The irrigation nozzle on the spindle should be aligned with the tool tip. In this occasion it should be noted that a uniform programming length is respected or the length of the 2 mm tool is estimated. The cooling medium air is sufficient.

5. Notes

■ Incorrect use of milling tools, polishing brushes, steam cleaner and water baths can overheat the material exposed and thus the properties are adversely affected. Deviations from the described manufacturing process can lead to different mechanical and optical properties of the material.
 ■ The lot number and the best before date are indicated on each Splint Memory packaging. In case of claims please always indicate the lot number of the product.
 ■ Do not use the product after expiration date.
 ■ The use of the material for the production of conventional splints is contraindicated, as is the production of denture bases. When properly handled and used by the dentist, orthodontist and dental technician, no side effects were observed. Undesirable biological reactions (e.g. allergies to material components) may occur in individual cases. Clarify the use of Splint Memory splints with your dentist or orthodontist in case of a known allergic response to components of Splint Memory.
 ■ Clean Splint Memory splints under cold running water with a toothbrush. Disinfect the splint with non-alcoholic detergents before first use.

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
Splint Memory milling blank for the generation of
occlusal splints with thermomemory effect

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6. Safety advice

■ Restorations are custom-made devices in accordance with the Directive 93/42/EEC. Record the BATCH no. for every procedure requiring identification of the material.
 ■ Storage: No special measures are required.
 ■ Disposal: Dental objects made of Splint Memory are not water soluble, are inert, and do not present a hazard for ground water. As such, they can be disposed of as domestic waste without any special precautions if not stated otherwise in your local disposal guidelines.
 ■ Warranty: Our recommendations concerning the application technique, regardless of whether they are communicated in writing, orally, or by means of practical instructions, are based on our own experiences and tests. As such, they are intended as guidelines only. We are continually striving to improve our products. Consequently, we reserve the right to make changes to their design and composition.
 ■ **Hazard Information:** Processing of Splint Memory blanks produces dusts which can irritate the eyes, skin, and airways. As such, it is essential to ensure the protective gear at your workplace is in perfect working order.
 ■ **ERNST HINRICHs Dental** is not liable for any damages caused by improper application of the material. To be used by trained specialist personnel for the purpose indicated only.

Manufacturer:
pro3dure medical GmbH
Am Burgberg 13 - 58642 Iserlohn, Germany
Phone: +49 (0)2374 920050-0 Fax: +49 (0)2374 920050-50

Distributor:
ERNST HINRICHs Dental GmbH
Borsigstr. 1 - DE-38644 Goslar
Phone: +49 (0) 53 21 - 37 90 - 0 Fax: +49 (0) 53 21 - 38 96 32
info@hinrichs-dental.de · www.hinrichs-dental.de
Made in Germany

Use by date
Haltbarkeitsdatum
Utiliser par date
Utilizar fecha de validez
Usare entro la data

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