# ERNST HINRICHS GmbH

Borsigstraße 1 · DE-38644 Goslar/Germany  
Telefon +49 (0) 53 21/5 06 24-25 · Telefax +49 (0) 53 21/5 08 81  
info@hinrichs-dental.de · www.hinrichs-dental.de

## Zirkon BioStar (weiß opak, Colour) Zirkon BioStar Z (weiß transluzent) Zirkon BioStar HT (hoch transluzent weiß, HT Colour) Zirkon BioStar S

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#### Name and general description

Porcelain frames type II, class 6 accord. DIN EN ISO 6872.

All **Zirkon BioStar** materials are dental materials (semi-finished products) made of yttrium stabilized, pre-sintered zirconium dioxide for milled production of crowns and bridge frameworks on commercial CAD/CAM systems or hand-operated copy-milling machines with outstanding biocompatibility and high resistance against tension and pressure.

The pre-sintered blanks are eminently suitable for machining and have excellent edge stability. On account of the special production process, extremely constant firing shrinkage is achieved, even over different batches, meaning that in most cases it is not necessary to adjust the enlargement factor on the CAD/CAM unit. Only diamond-charged grinding tools or hard metal tools are to be used for machining.

#### Material types

The material types differ regarding the strength value achievable after the final sintering and contain a different amount of aluminum oxide. They differ in the manner of whiteness and translucency and both show a different resistance against hydrothermal aging. **Zirkon BioStar Colour** (500, 800 1000, 1333, 2000) and **Zirkon BioStar HT Colour** (100, 200, 500, 800, 1000) are already persistent coloured in the pre-sintered stage. Consequently, the dental surgeon or dental technician can select the ideal base material for the respective case. On account of the high strength values, both materials are also suitable for multiple-element posterior bridgework.

**Zirkon BioStar HT** is a raw material with a reduced  $Al_2O_3$  content and a special initial grain size. Due to this arrangement this raw material produces an optimal relation of maximum translucency and an excellent hydrothermal consistency. Due to a higher density before sintering, this raw material is also suitable for the production of fully anatomical crowns, not to be layered over after the sintering process, but just slightly coloured before sintering instead.

If required, additional physical and chemical key values are referred in our material data sheets. All material types were tested by an accredited testing laboratory regarding biological compatibility and classified as eminently suitable for use in the dental sector.

#### Indication

All **Zirkon BioStar** products are suitable for the construction of dental restorations including crowns and bridges (both crowns and pontics as well as fully anatomic restorations with a maximum of two pontics in lateral areas and four pontics in anterior areas) or milled substructures, telescopic crowns, Maryland bridges, inlays, onlays, open end bridges (with a maximum of two pontics) and individual implant abutments.

#### Contraindications

In the case of bruxism, inadequate hard tooth structure or insufficient space, the use of full ceramic crowns and bridges is not recommended or only advised in some cases. The same applies in the case of inadequate oral hygiene or inadequate coordination of masticatory movements.

A full ceramic crown or bridge may also not be used in the event of unsatisfactory preparation results or in the case of production faults.

#### Preparation information

The minimum wall strength values in the following table are merely reference values. The decision on the manner and undertaking of the preparation is the responsibility of the dental surgeon and dental technician.

**The wall strength should be at least:**

Posterior teeth	occlusal 1.5 mm	circular 0.5 mm
Anterior teeth	incisal/occlusal 0.7 mm	circular 0.5 mm
Abutment crowns of bridge frameworks with 2 pontics	incisal/occlusal 1.0 mm	circular 0.7 mm

**Minimum connector areas:**

Anterior bridge framework	with one pontic	7 mm <sup>2</sup>
Anterior bridge framework	with two pontics	9 mm <sup>2</sup>
Posterior bridge framework	with one pontic	9 mm <sup>2</sup>
Posterior bridge framework	with two pontics	12 mm <sup>2</sup>

**Milling**

Milling can be undertaken with most commercially available CAD/CAM milling machines or hand-operated copy-milling machines. The machine manufacturer's manual should be consulted regarding the selection of the milling cutters and the necessary revolutions/minute. In order to avoid staining or a loss of strength, the blank should if possible not come into contact with liquids (water, oils) or fats (handcream) during milling and related handling of the blank.

Wet machining of the dental blanks in their pre-sintered state is in principle possible, but does require that the restoration be dried completely before sintering. Drying can be undertaken in a drying chamber at approx. 80°C or under an infrared lamp. Alternatively, a low-temperature firing program can be selected which will not damage the blank.

Notches, sharp edges or slight damage, also on the underside of the framework, may reduce the strength of the completed sintered work. Consequently, it is already important to ensure that surfaces are well rounded and smooth when milling.

**Colouring**

If required, the frameworks may be matched to the desired tooth colour with the standard colouring solutions which are licensed for zirconium oxide ceramics. In order to achieve even results, it must be ensured that the framework is completely clean, fat-free and dry before colouring. The instructions of the manufacturers of the colouring solutions should be followed for the further work procedures.

**Sintering**

The sintering of the preparations can be undertaken in all standard furnaces which are suitable for zirconium oxide ceramics. The manufacturer's operating instructions must be followed in all cases.

When sintered, the frameworks shrink by around 20 %. The required movement is only possible if the individual crowns or frameworks are fired on corresponding sinter bearings. The framework must be positioned in such a way that it cannot come into contact with other frameworks or the sinter carrier in order to avoid any deformation. The sinter carriers must be examined at least every two weeks in accordance with the manufacturer's instructions to ensure that they are level. Anterior and posterior dental frameworks with three or more elements should generally be positioned with the occlusal side upwards. If sintering wires or pegs are used, it must be ensured that there is no contact with the cap wall area. Alternatively, sintering may be carried out in a sintering crucible in which the restoration is embedded in suitable granular material. One suitable type of granules are ZrO<sub>2</sub> sintering spheres. When using sintering crucibles and granules, the framework is placed in the granules with the occlusal side downwards. The respective firing cycles can be completed in 8 to 12 hours. The dental laboratory must check and determine in the individual case whether these high heating rates are possible in the respective sintering furnace and whether they are still tolerated by the respective restoration without distortion or cracking.

The shrinkage of the preparations made of Zirkon BioStar and Zirkon BioStar Z remains at an extremely constant level, both within the individual batches and over several batches, on account of the special production process. We shall advise you of the exact values with each delivery.

**Please refer to the information/parameters in the instructions for use issued by the manufacturer of the equipment you are using.**

**Zirkon BioStar HT (high translucent)****Sintering:**

To achieve an optimum translucency, the use of a fresh sintering bowl as well as sintering ZrO<sub>2</sub>-beads is strongly recommended for the process of sintering Zirkon BioStar HT (high translucent). The sintering bowl must not have been used for sintering other zirconia materials previously.

**Colouring:**

Zirkon BioStar HT blanks may due to their higher density before sintering, be stained before the sintering process using appropriate stains and are therefore suitable for the production of fully anatomical and monolithic restorations without any over-layering or veneering.

## Zirkon BioStar S - Additional information for the instructions for use

### Instructions for the use of Zirkon BioStar S in your Sirona inLab® or inLab® MCXL system

#### General information:

Basically, you can use your Zirkon **BioStar S** blocks as other blocks for your Sirona® system.

#### Programming of the Sirona inLab® or inLab® MCXL System:

1. Place the Zirkon BioStar S block as usual.
2. Chose from the selection of blocks "**VITA In-Ceram®YZ**" the corresponding block size for the process.
3. Enter the "Z-Code" indicated on the package **completely**, including "\*\*\*", as for instance: **\*Z24912A**.
4. Inicie milling

VITA In-Ceram®YZ, Sirona inLab® and inLab® MCXL Systems are registered trademark of the manufacturers.

#### Please, consider when using BioStar HT Colour:

**Either sinter using an final temperature of 1.450 °C with a time of maintenance of 3.5 h, or 1.550 °C for 30 min, alternatively. All others values match the following chart. The achievable effects of translucency are not affected by either method.**

#### Depending on the material, the recommended sintering temperatures are:

Sintering program for Zirkon BioStar white opaque / HT (white) / Colour / S						
	single crown		bridge of 3 - 4 pieces		brigde of > 4 pieces & circular bridges	
	°C	min.	°C	min.	°C	min.
Preheating rate up to 900 °C °C / min.	15	59	12	73	8	110
Holding time at 900 °C	900	5	900	15	900	30
Preheating rate to final temperature of 1.450 °C	15	36	12	46	8	69
Holding time at final temperature of 1.450 °C	1450	30	1450	120	1450	120
Cooling rate down to 900 °C °C / min.	10	55	5	110	4	138
Holding time at 900 °C	900	5	900	10	900	20
Cooling rate down to 200 °C	10	70	5	135	4	175
Total time in h	4 h 20 min.		8 h 29 min.		11 h 02 min.	

Sintering program for Zirkon BioStar Z white translucent						
	single crown		bridge of 3 - 4 pieces		brigde of > 4 pieces & circular bridges	
	°C	min.	°C	min.	°C	min.
Preheating rate up to 900 °C °C / Min.	15	59	12	73	8	110
Holding time at 900 °C	900	5	900	15	900	30
Preheating rate to final temperature of 1.500 °C	15	40	12	50	8	75
Holding time at final temperature of 1.500 °C	1500	30	1500	120	1500	120
Cooling rate down to 900 °C °C / min.	10	60	5	120	4	150
Holding time at 900 °C	900	5	900	10	900	20
Cooling rate down to 200 °C	10	70	5	135	4	175
Gesamtzeit h	4 h 29 min.		8 h 43 min.		11 h 20 min.	

As the temperature measurement in every furnace can give somewhat different results, it may possibly be necessary in an individual case to adjust the specific furnace temperature by carrying out a test firing cycle. Massive drops in temperature when opening the sintering furnace should be avoided, as opening the heated furnace can have a negative effect on both the furnace and the sintered products, damaging the restoration. The temperatures vary according to furnace type. With our zirconia products, avoid opening the furnace at temperatures above 250°C. Please refer to the instructions provided by the furnace manufacturer. Dental frameworks and crowns made of zirconium oxide should be fired in a furnace which is only used for these products. If additionally restorations made of other ceramic materials are fired or if liquid ceramics are baked in the same furnace, this may impede the sintering behavior or result in localized discolouration. If the firing of different materials in the same furnace cannot be avoided, then an empty firing cycle at maximum temperature is recommended before firing zirconium oxide ceramics.

### **Facing of the frameworks**

Crowns and bridge frameworks made of Zirkon BioStar materials may be faced with standard facing ceramics for frameworks made of zirconium oxide. The facings passed by the respective manufacturers for zirconium oxide are adapted regarding the thermal expansion coefficient and the material bonding. Please consult the descriptions of the manufacturers of the facing material for information regarding whether thermal or chemical pretreatment of the sintered construction is recommended before facing.

### **Finishing after final sintering**

Any necessary finishing of the sintered constructions must only be done using water cooling during grinding operations, thus avoiding localized overheating and damage to the material.

### **Information on storage**

Store the zirconium discs at normal room conditions and humidity ("cool and dry" corresponding 10 °C – 35 °C and 30 up to 70 % relative humidity). Even storing the material for a longer period, we do not expect a change of the technical characteristics. Avoid a direct insolation.

**The minimum wall thicknesses for crown and bridge frameworks and sintering program heating rates given in the instructions for use have been tested many times in practice and have proven satisfactory. However, it may be necessary due to unusual framework size or unequal material concentrations to adjust the processing conditions to the specific restoration. Therefore simply adhering to the listed parameters cannot guarantee accurate results.**

### **Safety hint**

The inhalation of ceramic dust is harmful to your health. Avoid contact with skin and eyes. When preparing and blasting use a suction extraction system and breathing mask!

### **Disposal of material residues**

Remains of the blanks and, of course, also sintered parts made of Zirkon BioStar materials are categorized as chemically completely inert and are not subject to any special disposal stipulations.

### **Secondary effects**

Such as allergies to zirconium dioxide may very rarely occur.

### **Reciprocal actions**

unknown.

### **Reactions**

In case of known incompatibilities and allergies to contents of the zirconium dioxide.

**Distributor:** ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar

**Manufacturer:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar

SILADENT Dr. Böhme & Schöps GmbH carries the mark  0124.

## Material data sheets for the different Zirkon BioStar materials:

<b>Chemical formula:</b>	ZrO <sub>2</sub> / Y <sub>2</sub> O <sub>3</sub> / Al <sub>2</sub> O <sub>3</sub>
<b>Chemical name</b>	Ytria-stabilized tetragonal zirconia
<b>Colour:</b>	white opaque / white translucent / tooth coloured in various shades / precoloured high translucent
<b>Product description:</b>	Presintered zirconium oxide components for dental use

<b>Chemical composition</b>		Zirkon BioStar Z	Zirkon BioStar	Zirkon BioStar Colour	Zirkon BioStar HT white	Zirkon BioStar HT Colour
ZrO <sub>2</sub> / HfO <sub>2</sub> / Y <sub>2</sub> O <sub>3</sub>	wt%	> 99,0	> 99,0	> 99,0	> 99,0	> 99,0
Y <sub>2</sub> O <sub>3</sub>	wt%	5,15 +/- 0,20	5,15 +/- 0,20	5,15 +/- 0,20	5,15 +/- 0,20	5,65 +/- 0,30
HfO <sub>2</sub>	wt%	< 5,0	< 5,0	< 5,0	< 5,0	< 5,0
Al <sub>2</sub> O <sub>3</sub>	wt%	< 0,1	< 0,3	< 0,3	< 0,1	< 0,1
Fe <sub>2</sub> O <sub>3</sub>	wt%	< 0,1	< 0,1	< 0,3	< 0,1	< 0,1
Na <sub>2</sub> O	wt%	< 0,04	< 0,04	< 0,04	< 0,04	< 0,04

<b>Physical material data after final sintering</b>						
Density	g/cm <sup>3</sup>	> 6,05	> 6,05	> 6,00	> 6,05	> 6,05
Bending strength	MPa	1.000 +/- 200	1.200 +/- 200	1.000 +/- 200	1.000 +/- 200	1.000 +/- 200
Modulus of elasticity	Gpa	210	210	210	210	210
Thermal expansion	10 <sup>-6</sup> 1/K	10,7	10,7	10,7	10,7	10,7
Open porosity		none	none	none	none	none

<b>Radioactivity</b>						
Activity <sup>238</sup> U, <sup>226</sup> Ra, <sup>232</sup> Th	Bq/g	< 0,2	< 0,2	< 0,2	< 0,2	< 0,2

*meets the requirements of ISO 6872 and ISO 13356.2*

### **Biocompatibility and chemical solubility:**

Cytotoxicity	Non-toxic	DIN EN ISO 10993 / DIN EN ISO 9363
Chemical solubility	below detection limit (< 1 µg/cm <sup>2</sup> )	DIN EN ISO 10903 / DIN EN ISO 6872

The above data are typical values and relate only to the test pieces from which they were obtained. They are therefore not unconditionally applicable to other forms of the same material. The values are therefore to be taken as indicative only.

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**Commerciale:** ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
**Fabbricante:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH segnare il confezionamento a ☎ 0124.