Hinri-Alloy CB

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Hinri-Alloy CB

Instructions for use

Area of application: Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 5.

Composition (in % by mass):

Co: 64 Cr: 21.0 Mo: 6.0 W: 6.0

Other constituents: Si, Fe, Mn

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines):

Proof stress 0.2 % (MPa):	570
Elongation at rupture (%):	10
Modulus of elasticity (GPa):	194
Coefficient of expansion:	20°C-500°C 14.1 x 10 ⁻⁶ /K 20°C-600°C 14.6 x 10 ⁻⁶ /K
Tensile strength (MPa):	734
Vickers hardness VH 10:	286
Density (g/cm ³):	8.8
Melting range (°C):	1.309 - 1.417
Casting temperature (°C):	1.460
Max. oxide firing temperature (°C):	935

Product description: Hinri-Alloy CB is a non-precious, cobalt based bonding alloy. Cause of the technical properties very thin-walled constructions are possible. The structure of the alloy generates less oxidation also after repeated casting and offers a gut processing and polishing. Its thermal expansion coefficient is ideal for every ceramic of the latest generation. Hinri-Alloy CB is highly corrosion resistant and does not contain any beryllium, indium or gallium ($\leq 0.1\%$).

Investing and casting: Hinri-Alloy CB can be used with phosphate-bonded crown and bridge investments, such as Hinrivest KB. It is preheated to 850 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate unglazed and preheated ceramic crucible for casting Hinri-Alloy CB to prevent it being contaminated by other alloys. Clean the crucible after each cast. With the high-frequency melting the alloy is cast when all the cylinders have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. Please note that the melting point may vary depending on the type of casting machine used and when melting with vacuum, for example, less oxide forms and the oxide layer in the vacuum disintegrates more quickly. In general cast with incorporating flux. Considerable experience is required when melting high-grade alloys and it is most important that heating is set exactly according to the manufacturer's instructions. Sample values for flame regulation: Acetylene 0,4 bar / Oxygen 2 bar; Propane 0,2 bar / Oxygen 2 bar; Leading pressure methane / Oxygen 2 bar. Inaccurate heat setting impairs the alloy and can have a detrimental effect on the porcelain bond. Place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multiflame welding torches only. Do not use any flux. Casting Hinri-Alloy CB repeatedly is not recommended, as this could alter both its composition and properties. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

Manufacturing and polishing: Cool down the muffle at room temperature, remove the investment and sandblast with aluminum oxide 100 µm to 250 µm at a pressure of 3-4 bar. Continue manufacturing process using fraises for hard metal.

Firing porcelain:

1. The worked framework has to be sandblasted with disposable equipment in aluminum oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steamclean. After cleaning, the framework should not be touched anymore by hand.
2. Fire for 5 min. under vacuum at 950 – 980 °C (10° C more than back-ground material firing temperature).
3. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminum equipment 110 to 150 µm at a pressure of 2.5 - 3 bar, then steam.
4. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.
5. It is possible to use all ceramics normally available on the market for CoCr alloys. Apply manufacturing instructions of the ceramic producer (e.G. VITA VM13). Do not use any bonder.

Soldering and welding: Presolder using ERNST HINRICHS CrCoMo special solder (REF 893060 / 893061). For laser welding use ERNST HINRICHS CrCo laser welding rod (REF 893062).

Safety hints: Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask. The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. A Patch-Test is thus advisable. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

Distributor: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar

Manufacturer: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia

Mesa Italia Srl. packaging carries the mark

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Commerciale: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travaagliato / Brescia
Mesa Italia Srl. segnare il confezionamento a
ERNST HINRICHS Dental GmbH