Hinri-Alloy CoCr

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Hinri-Alloy CoCr

Instructions for use

Area of application: Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5

Product description: Hinri-Alloy CoCr is an alloy for denture bases is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame (0.5 mm in wax; 0.3 mm in metal). Moreover its molecular structure allows to obtain smooth, compact surfaces with little oxide formation. It is highly corrosion resistant and is free of toxic elements like beryllium, indium and gallium ($\leq 0.1\%$).

Composition (in % by mass):

Co: 62.0 Cr: 31.0 Mo 5.0

Other constituents: Si, C, Mn

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines):

Proof stress 0,2 % (MPa):	650
Elongation at rupture (%)	5.0
Modulus of elasticity (GPa):	220
Tensile strength (MPa):	890
Vickers hardness VH 10:	350
Density (g/cm ³):	8.3
Melting range (°C):	1.280 - 1.360
Casting temperature (°C):	1.450

Investing and casting: Hinri-Alloy CoCr can be used with phosphate-bonded CrCo investments such as Hinriinvest® CoCr or Hinriinvest® SG. Prepare the moulds for casting according to the ERNST HINRICHS ring less system. Preheating range is 850 °C – 1.000 °C. ERNST HINRICHS recommend 920 °C. Heat soak at casting temperature for 30 minutes.

The alloy is cast with an induction machine when all the cylindrical ingots have melted together and a final shadow runs over the molten metal. Considerable experience is required when melting high-grade alloys automatically, and it is most important that the heating is set exactly. Accurate heat setting is essential to protect the properties of the alloy. Light the flame and adjust the pressure in order to have a blue nucleus about 4 mm long. Place some metal ingots into the preheated crucible. Move the flame over the crucible rotating the tip, in order to evenly heat the alloy. As soon as the ingots are wholly melted, and after disappearing of their shadows, start centrifugal action. Casting Hinri-Alloy CoCr repeatedly is not recommended. During re-use pay attention that the alloys are from the same batch. Mix old and new alloys in a ratio of 1:1.

Preparing and polishing: Cool down the muffle to room temperature, Sandblast the framework with aluminum oxide 120-150 µm and continue the working cycle according to your personal method. Remove carefully any small casting blows, go over all the corners and edges with a round bur and then smooth the surfaces with rubber abrasives and polish them.

Soldering and welding: Presolder using ERNST HINRICHS CrCoMo solder (2 mm, REF 893060, 1 mm, REF 893061). For Laser welding use ERNST HINRICHS CrCo laser (REF 893062) welding rod.

Safety hints: Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask. The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. A Patch-Test is thus advisable. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

Distributor: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar

Manufacturer: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia

Mesa Italia Srl. packaging carries the mark

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Commerciale: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia
Mesa Italia Srl. segnare il confezionamento a
ERNST HINRICHS Dental GmbH