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info@hinrichs-dental.de  
www.hinrichs-dental.de

## Hinri-Alloy N Instructions for use

**Area of application:** Non precious casting alloy based on nickel for metal ceramic accord. DIN EN ISO 22674, type 3.

**Product description:** Hinri-Alloy N is a non-precious, nickel based bonding alloy. Cause of the technical properties very thin walled constructions until 0.2 mm are possible. The structure of the alloy generates less oxidation also after repeated casting and offers a gut processing and polishing. Its thermal expansion coefficient is ideal for every ceramic of the latest generation. Hinri-Alloy N is highly corrosion resistant and does not contain any beryllium, indium or gallium ( $\leq 0.1\%$ ).

 **This alloy contains nickel and should not be processed by individuals with known nickel sensitivity (allergy). Metal dust containing nickel should not be inhaled. Hinri-Alloy N should not be used with anyone allergic to nickel.**

### Composition (in % by mass):

Ni: 63.0 Cr: 25.0 Mo: 9.0 Si: 2.0 Nb: 1.0

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

### Technical data (guidelines):

Proof stress (MPa):	360
Elongation at rupture (%):	26
Modulus of elasticity (GPa):	191
Coefficient of expansion:	25-500 °C $13.7 \times 10^{-6} \text{ K}^{-1}$ 25-600 °C $14.0 \times 10^{-6} \text{ K}^{-1}$
Tensile strength (MPa):	602
Vickers hardness VH 10:	180
Density (g/cm <sup>3</sup> ):	8.3
Melting range (°C):	1.298—1.344
Max. oxide firing temperature (°C)	950
Casting temperature (°C)	1.400

**Investing and casting:** Hinri-Alloy N can be used with phosphate-bonded crown and bridge investments, such as Hinrivest KB. It is preheated to 850 - 950 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate ceramic crucible for casting Hinri-Alloy N high-grade, non-precious alloy to prevent it being contaminated by other alloys. To avert shrinkage blisters and cracks avoid hollows, sudden elbows and impurities. Wax walls should be 0.5 mm thick, so that after manufacturing process metal wall thickness is at least 0.3 mm. Transitions from wax to metal should be realized with gradual passages. High-frequency melting equipments: When last ingots have melted in the preheated ceramic crucible and incandescence shadow has disappeared, start melting process immediately. Do not use any flux. Open-flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux. Sample values for flame regulation: acetylene 0.4 bar / oxygen 2 bar ; propane 0.2 bar / oxygen 2 bar; leading pressure methane / oxygen 2 bar. Overheating the melted material could cause shrinking hollows, microporosity and building of a coarse-grained structure, and cause bridges to break or ceramic investments to crack. Casting Hinri-Alloy N repeatedly is not recommended, as this could alter both its composition and properties. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

**Manufacturing and polishing:** Cool down the muffle at room temperature, remove the investment and sandblast with aluminum oxide 100 µm to 250 µm at a pressure of 3-4 bar. Continue manufacturing process using fraises for hard metal.

### Firing porcelain:

1. The surface has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steamclean. After cleaning, the framework should not be touched anymore.
2. Fire for 5 min. under vacuum at 950°-980°C (10°C more than back-ground material firing temperature).
3. After firing, the oxide layer has to be carefully sandblasted using disposable aluminium oxide equipment 110 to 150 µm at a pressure of 2.5 - 3 bar, then steam-clean.
4. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.
5. It is possible to use all ceramics normally available on the market for NiCr alloys. Apply manufacturing instructions of the ceramic producer (e.G. VITA VM13, Shofu Vintage Halo).

**Soldering and welding:** Presolder using ERNST HINRICHS CoNiCr solder (REF 893064) and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys. Laser weld using ERNST HINRICHS Cr-Co laser welding rod (REF 893062).

**Safety hints:** Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask! The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

**Contraindications, side effects and interactions with other dental alloys:** In cases of hypersensitivity (allergy) to the constituents of the alloy (a patch-test is thus advisable). Hypersensitivity reactions (allergies) or electrochemically-induced local dysaesthesia have been reported in individual cases. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

**Storage conditions:** No special storage conditions are required.

**Further questions:** Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

**Distributor:** ERNST HINRICHS Dental GmbH, Borsigstr. 1, DE 38644 Goslar

**Manufacturer:** Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia

Mesa Italia Srl. packaging carries the mark

CE 0123

Date of information: 09/2020

Printing Date: 16.09.2020

**ERNST HINRICHS Dental GmbH**

Borsigstr. 1  
D-38644 Goslar

Tel.: +49 (0) 53 21 - 5 06 24 / 25  
Fax: +49 (0) 53 21 - 5 08 81

info@hinrichs-dental.de  
www.hinrichs-dental.de