## Hinrisil® Speed

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#### Typ:

Additionsvernetzendes  
1 : 1 Dubliersilikon  
DIN EN ISO 14356,  
Typ 2 - irreversible  
Dubliermasse

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#### Technical Data:

Mixing ratio A (white) : B (yellow):	1 : 1
Mixing under vacuum:	40 sec.
Working time at 23°C:	> 3-4 min.
Curing time at 23°C:	10 min.
Tensile strength:	approx. 2.2 MPa
Elongation at break:	approx. 310%
Tear strength:	> 6.5 N/mm
Shore A hardness after 10 minutes:	24

#### Application:

Rapid-curing, addition-curing duplicating silicone, developed especially for express jobs using adhesive crepe sleeves or flasks. Hinrisil® Speed can be removed from the duplicating form after just 10 minutes and is suitable for use with the flaskless duplicating system. Models to be duplicated should be at room temperature and should not be soaked. Any dirt on the surface is removed by spraying with Hinrisil (REF 109400) and by blow-drying carefully. Materials, such as gypsum hardener, grease, oils and some separators, can produce a negative reaction with Hinrisil® Speed (refer to „Incompatible materials“ below). After blocking out with a light-curing composite such as Visio®-Bloc (ESPE), the greasy layer on the surface should be carefully removed using ethyl alcohol. Duplicating tape (REF 891007) is then placed round the model and a suitable stabilising plate selected. The retention holes in the plate should almost cover the dental arch and should not come into contact at all with the duplicating tape. The stabilising plate is secured in the fixation device (REF 891001). Mix Hinrisil® Speed under vacuum and pour it from a height of approx. 40 cm onto the centre of the model. Do not pour directly onto the contours of the model. Pour until the top of the model is completely covered. Set the model under the fixation device, place the stabilising plate in position and secure it approx. 5 mm above the dental arch. The fixation device should not be moved again until the silicone is cured. At 23 °C working temperature, the plaster model can be removed with the help of compressed air from the duplicating mould after 10 minutes. Carefully remove any flash or excess material with a pair of scissors. Lower working temperatures may lengthen the curing time by approx. 1 minutes; higher working temperatures may shorten the curing time respectively. Hinrisil® Speed can be measured using a dispenser because both components have almost the same viscosity.

#### Chemical characterisation:

Polymerised dimethyl polysiloxane filler mixture

#### Dangerous reactions:

No known dangerous reactions, if stored and used properly.

#### Toxicity:

None known

#### Directions:

Store containers tightly sealed in a well-ventilated place. Keep containers dry. Keep away from water, alkalis, strong acids and oxidising agents.

#### Ecological information:

Present experience suggests that there is no detrimental effect on the environment, if it is used properly and current regulations are observed. It is non-biodegradable and insoluble in water when cured. It can be easily separated from water by filtering.

#### Please note:

Cured silicone is not harmful to the environment and can be disposed of along with domestic waste, as long as this does not contravene local bye-laws.

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. ERNST HINRICHS products are subjected to strict quality controls.

We reserve the right to make technical changes.

#### Incompatible materials:

The following materials can produce a negative reaction with addition-curing silicones:

soft polyvinyl chloride, polysulphide MIL-S-8516, Mystik® 6207 adhesive tape, Mystik® 6215 adhesive tape, Scotch® cellophane strips, Scotch® 360 adhesive tape, Permacel® masking tape, Pliobond® adhesive, Neoprene rubber Buna N rubber, Vitron® A rubber, colophony-based soldering flux, acid-based soldering flux, disinfectant, ISOLIT, cyanoacrylate glue, amino-hardened epoxy, Humiseal® 6215 adhesive tape, vinyl insulating tape, latex vacuum hoses, GRS rubber, natural gum, certain types of silicone and RTV rubber, SUPER-SEP, milling oils, sulphur compounds such as thiols, sulphides, sulphates, sulphites, thiourea and nitrogen compounds such as amines, amides, imides, acids.

If there is any doubt, we recommend testing the material before use.

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