HINRIVEST® uni 3D

DIN EN ISO 15912, Typ 1, Klasse 1 & 2

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Instructions for use

Description:

HINRIVEST uni 3D is a phosphata-bonded, graphite-free speed casting investment with exceptional properties. It is suitable for casting all dental alloys. It can be used for crowns and bridges made of precious and non-precious alloys and press ceramics (lithium disilicates). Furthermore, crowns and bridges printed in 3D resin or milled from PMMA can be cast.

Technical data:

Mixing liquid:	Hinrivest® type 100 expansion liquid
Mixing ratio, Powder : liquid:	100 g : 25 ml
Mixing under vacuum:	60 sec.
Working time:	5 - 6 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	9 - 11 min.

Modelling: The best results in terms of fit and surface finish are achieved with pure wax models. Resins can lead to poorer surfaces and muffle cracks during high-speed casting.

Mixing: Mix the investment material by hand with a spatula until an even coating is achieved. Place under vacuum for 15 seconds without switching on the agitator, then stir for 60 seconds. Check the vacuum agitator constantly for function. An insufficient vacuum results in fitting tolerances and bubbles on the cast object.

Processing time: The processing time is approx. 5 to 6 min. at approx. 22°C (room temperature). Caution: The processing time depends on the room temperature. Heat shortens the processing time.

Investing: Filling the muffle with investment material: The vibrator should only be used to assist the flow if it appears necessary. Avoid vigorous vibration! This can lead to the formation of bubbles and segregation of the investment material.

Placing the muffle in the furnace: After 20 minutes, calculated from the start of mixing! Roughen the top of the muffle before placing it in the furnace!

Attaching the casting sprues: The cast objects are placed in the muffle in such a way that they are never in the thermal centre. The thin parts (e.g. crown margins) point towards the muffle wall, distance to the muffle edge 5 mm.

Concentration of the expansion liquid:

Liquid concentration (1.000 ml liquid)			Concentrations (Liquid / Water)	
Concentration	Concentrate (ml)	dist. water (ml)	100 g	160 g
50%	500	500	12.5 ml / 12.5 ml	20 ml / 20 ml
60%	600	400	15 ml / 10 ml	24 ml / 16 ml
65%	650	350	16 ml / 9 ml	26 ml / 14 ml
70%	700	300	17.5 ml / 6.5 ml	28 ml / 12 ml
75%	750	250	19 ml / 6 ml	30 ml / 10 ml
95%	950	50	22.5 ml / 2.5 ml	36 ml / 4 ml
100%	1.000	-	25 ml / -	40 ml / -

	Liquid ratio / distilled water:	100 g	160 g
High-gold-content precious metals	60 % / 40 %	15 ml / 10 ml	24 ml / 16 ml
Semi-precious metals	65 % / 35 %	16 ml / 9 ml	26 ml / 14 ml
Bonding alloys	70 % / 30 %	18 ml / 7 ml	28 ml / 12 ml
Semi-precious bonding alloys	75 % / 25 %	19 ml / 6 ml	30 ml / 10 ml
Precious free alloys (CrCo)	90 % / 10 %	23 ml / 2 ml	36 ml / 4 ml

Speed casting: HINRIVEST® uni 3D is preferably preheated in the speed casting method. Preheat the furnace to 850 °C and, exactly 20 minutes after mixing the investment, place the mould in the furnace with the sprue hole facing downwards. (Adhere strictly to the prescribed time, which is measured from the first contact of the powder and liquid.) Moulds should not be placed directly onto the floor of the furnace. We recommend using either a grooved tray for collecting wax.

Caution! during the first 15 minutes the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. If a preheating temperature other than 850°C is required, the temperature can be adjusted 15 minutes after placing the mould in the furnace. Moulds are ready to cast after heat soaking for 45 - 60 minutes depending on the size of the ring.

Conventional preheating: If a conventional preheating is necessary should allow the mould to set for at least 60 minutes, then place it in a cold furnace and heat to the required temperature at a heat rate of 5 - 7°C per minute. Holding stages are not required when using HINRIVEST® uni 3D.

Press ceramic: For the production of press ceramic objects, a liquid concentration (expansion liquid type 100) of 65% is recommended. Due to the good controllability of the investment material, the expansion (higher or lower expansion) can be adjusted or controlled depending on the situation.

Please note: Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. ERNST HINRICHS Dental products are subjected to strict quality controls. We reserve the right to make technical changes.

Warning: Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

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