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
## Keralloy® BioStar

### Instructions for use

**Area of application:** Non precious milling alloy based on cobalt for metal ceramic accord. DIN EN ISO 22674, type 4.

**Product description:** Keralloy® BioStar is a non-precious, cobalt based bonding alloy. Cause of the technical properties very thin walled constructions until 0.2 mm are possible. The structure of the alloy generates less oxidation also after repeated casting and offers a gut processing and polishing. Its thermal expansion coefficient is ideal for every ceramic of the latest generation. Keralloy® BioStar is highly corrosion resistant and does not contain any beryllium, indium or gallium ( $\leq 0.1\%$ ).

#### Composition (in % by mass):

Co: 60.0 Cr: 28.0 W: 9.0 Si: 1.5

**Other constituents:** Mn, Fe

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

#### Technical data (guidelines):

Proof stress (MPa):	361
Elongation at rupture (%):	16
Modulus of elasticity (GPa):	183
Coefficient of expansion:	25-500 °C $14.2 \times 10^{-6} \text{ K}^{-1}$ 25-600 °C $14.4 \times 10^{-6} \text{ K}^{-1}$
Tensile strength (MPa):	456
Vickers hardness VH 10:	273
Density (g/cm <sup>3</sup> ):	8.5
Melting range (°C):	1.308 -1.386
Max. oxide firing temperature (°C)	980
Casting temperature (°C)	1.436

**Detach of the Framework:** Detach the framework and single parts with suitable hard metal mills or cutting wheels, and slur/grind.

**Manufacturing and polishing:** Cool down the muffle at room temperature, remove the investment and sandblast with aluminum oxide 100  $\mu\text{m}$  to 250  $\mu\text{m}$  at a pressure of 3-4 bar. Continue manufacturing process using fraises for hard metal.

#### Firing porcelain:

1. The surface has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150  $\mu\text{m}$  at a max. pressure of 3-4 bar and then steamclean. After cleaning, the framework should not be touched anymore.
2. Fire for 5 min. under vacuum at 950°-980°C (10°C more than back-ground material firing temperature).
3. After firing, the oxide layer has to be carefully sandblasted using disposable aluminium oxide equipment 110 to 150  $\mu\text{m}$  at a pressure of 2.5 - 3 bar, then steam-clean.
4. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.
5. It is possible to use all ceramics normally available on the market for NiCr alloys. Apply manufacturing instructions of the ceramic producer (e.G. VITA - VMK Master®, Elephant Dental - Interaction Antagon). Do not use any bonder.

**Soldering and welding:** Presolder using CrCo solder and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys.

**Safety hints:** Metal dust and smokes are dangerous to your health. When deflasking and blasting use a suction extraction system and breathing mask! It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

**Contraindications, side effects and interactions with other dental alloys:** In cases of hypersensitivity (allergy) to the constituents of the alloy (a patch-test is thus advisable). Hypersensitivity reactions (allergies) or electrochemically-induced local dysaesthesia have been reported in individual cases. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

**Storage conditions:** No special storage conditions are required.

**Further questions:** Please contact our technical services (Tel.: +49 5321-50624 / 25) or our sales representatives.

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**Commerciale:** ERNST HINRICHS Dental GmbH, Borsigstr. 1, DE 38644 Goslar
**Fabbricante:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH segnare il confezionamento a **CE 0124**
**ERNST HINRICHS Dental GmbH**