PMMA BioStar Gebrauchsanweisung

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PMMA BioStar

Instructions for use

Product description: Dental milling discs based on PMMA (polymethyl methacrylate) which burn out residue-free and are developed for the casting technique. PMMA BioStar is available in 3 different colours (transparent, blue, ivory) and in the sizes of 14, 18, 25 and 30 mm (Ø 98.5 mm with shoulder).

Technical Data:

Chemical characterization	100 % PMMA (polymethyl methacrylate)
Density	1.18 g/cm ³
Ball indentation hardness	175 MPa
E-Modulus	3.300 MPa

General information: Please note that the improper use of milling tools, polishing brushes, steam jets and water baths can cause the material to overheat. Do not exceed temperatures of 150 °C. This can result damage to the material.

Milling parameters: Follow the instructions of the manufacturers of CAD-CAM machines and strategies. For the dry milling of PMMA BioStar please use carbide metal milling tools. To cut out the milled restorations from the PMMA BioStar use burs or discs for cutting acrylic materials. Avoid temperatures of more than 150 °C, which can damage the disc.

Cleaning: If a cleaning process is necessary, clean it manual or use an ultrasonic cleaner with a short setting of about one minute at a temperature of no more than 40 °C. Cleaning solutions should be added only in low concentrations.

Attached sprues: Please check the rules for your casting alloys or pressable ceramics to design sprues and runners for casting. Acrylic molds for casting should be covered with a thin layer of wax before embedding.

Embedding: Use a ring free embedding system (e. G. ERNST HINRICHS system) or coat the ring with a ceramic ring liner before filling with investment. Mix the investment material using vacuum and according to instructions for use of the manufacturer. Heat the mold conventionally or with a speed heating system up to 580 °C, holding time 30 min. Afterwards heat the mould to the required end temperature according to the used alloy or press ceramic. Follow the manufacturer instructions for use of the investment material, the alloy and / or press ceramic used.

Disposal of waste: The PMMA BioStar product is made from pure PMMA. Waste can be recycled with other acrylics or can be disposed of through domestic residual waste system. PMMA is insoluble in water and is no hazard to groundwater.

Further questions: Please contact ERNST HINRICHS technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

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