CoCr BioStar

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
CoCr BioStar

Instructions for use

Area of application: Non precious milling alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 4.

Product description: CoCr BioStar is non precious milling alloy based on cobalt for metal-ceramic used in dental milling machines (CAD-CAM). It does not contain any nickel, beryllium or gallium. One of the remarkable features is the high corrosion resistance and biocompatibility. CoCr BioStar is suitable for soldering.

Composition (in % by mass)

Co	61.65
Cr	27.75
W	8.45
Si	1.61

Other constituents:

Mn-Fe

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines)

Proof stress 0.2 % (MPa)	375
Density g/cm ³	8.55
Vickers hardness VH 10	285
Coefficient of expansion 20°C – 600°C	14.5 x 10 ⁻⁶ K ⁻¹
Melting range °C	1.400 – 1.418
Modulus of elasticity GPa	240
Elongation at rupture A5 %	14.7
Casting temperature (°C):	1.518
Max. oxide firing temperature (°C):	980

Detach of the Framework: Detach the framework and single parts with suitable hard metal mills or cutting wheels, and slur/grind.

Firing porcelain:

1. An oxide firing (atmospheric at 980 °C, 10 minutes) is not absolutely mandatory but could have advantages.
2. After the oxide firing blast always the oxide layer with fresh aluminum oxides 110 µm and clean the surface with a steam-jet. Do not clean with hydrogen fluoride!
3. The use of a bonder is positive for the metal-ceramic bonding. Wash and opaque firing according the manufacturer's instructions for use.
4. We recommend a slow cooling depending to the size of the construction according the manufacturer's instructions for use (Reference: Vita VM 13 / VMK Master Bonder).

Cleaning: Clean frames with steam or ultrasonic using distilled water. **Caution:** Do not strip chemically.

Soldering and welding: Presolder using CrCo solder and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys.

Polishing: Remove residues of oxides with glass beads. Glaze, gum and polish the elaborated metal surfaces with suitable abrasive tools.

Recommendation for single use: A recycling of the alloy is not recommended.

Safety hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask type FFP2!

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Storage conditions: No special storage conditions are required.

Further questions: Please contact ERNST HINRICHS technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

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Distributor: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Manufacturer: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH packaging carries the mark
ERNST HINRICHS GmbH