CoCrMo Spezial-Lot

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CoCrMo Special Solder

Instructions for use

Area of application: Non precious cobalt based hard solder alloy acc. EN ISO 9333

Description: CrCoMo special solder is indicated to weld cobalt based alloys. It's investible with ceramic and completely without nickel, beryllium, gallium and indium ($\leq 0.1\%$).

Working instruction: Mix a small amount of the added flux powder with a drop of water and cover the warmed up soldered part. Gap must not exceed 0,2 mm, heat the area to be soldered, apply flux, add solder and heat by a flame. In some cases you can add a small amount of solder drops directly on the position which has to be soldered. The highmelting point guarantees that the solder has a high diffusion depth. It can be used for non visible enlargements and reparations. If you handle the special solder correctly the quality will be the same as welding our laser process.

Composition (in % by mass): Co: 62.0 Cr: 29.0 Mo: 4.0 Si: 3.0 Other constituents: C, Mn, Fe

Melting range (°C): Solidus-liquidus temperature: 1.071 - 1.260 °C Flow temperature: 1.320 °C

Safety hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask!

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy. Hypersensitivity reactions (allergies) or electrochemically-induced local dysaesthesia have been reported in individual cases. A Patch-Test is thus advisable. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. The alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24) or our sales representatives.

Distributor: ERNST HINRICHS Dental GmbH, Borsigstr. 1, DE 38644 Goslar

Manufacturer: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia

Mesa Italia Srl. packaging carries the mark

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El fabricante: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travaigliato / Brescia
Mesa Italia Srl. rotula el envase comercial con
Mercante: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: Mesa Italia Srl., Via dell' Artigianato, 35/37/39, IT 25039 Travaigliato / Brescia
Mesa Italia Srl. contrassegna la confezione con il marchio