HinriDon

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HinriDon

Instructions for use

HinriDon is an economic, cadmium-free denture base resin that guarantees high quality results by using any known hot-curing polymerization technique. HinriDon belongs to the newest generation of short-term hot-curing dental acrylics and also perfectly suitable for the injection technique.

Technical data:

Mixing ratio Powder/Liquid:	10 : 4
Thickening phase approx.	10 min.*
Processing- and pressing phase approx.	45 min.*
Polymerization at 95 °C	20 min.*
Flashpoint:	10 °C
Max. residual monomer content:	0,8 %

* The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures extend, lower temperatures shorten these processing times.

Indications:

- Total upper and lower dentures by using the pressing-/packing technique

Product features: Cadmium-free, high mechanical properties, physiologic shades and 100% color stable, high quality results and an accurate fitting denture guaranteed, easy handling, mucous membrane friendly due to a low amount of residual monomer, compatible to all known individual shades, evaluated and certified biocompatibility. HinriDon allows quick and easy handling with no fancy equipment needed.

HinriDon does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded), supplied with retentions and wetted a bonding agent.

Preparation: The wax modellation is embedded into the gypsum (e. g. Hinridur). The flask is boiled out thoroughly to remove the wax before applying alginate based isolation on the surface of the gypsum. Hinrisep K always shows reliable results. **Important:** The isolation needs to be well-dried to make sure there are no white marks or spots on the finished acrylic.

Processing: Place required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as needed to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a swelling phase of approximately 5 to 10 min. the material can be handled within another 45 min.

Pressing: Place the material with a little excess into the flask. Slowly close the flask and continuously apply pressure. Keep the pressure for at least 10 min. A test pressing is not necessary due to the excellent flow ability of the acrylic.

Suggested polymerization: This acrylic is especially designed for the short-term hot-curing polymerization. Immediately place the flask with a clamp in boiling water (approx. 3 l of water per flask) and allow curing for about 20 min. For curing larger projects place the flask in boiling water, turn heat source off and wait for 15 min. Bring water to a boil again and allow curing for 10 min. Of course the material can be processed by using any other polymerization technique as well.

After completion: Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste (0,8% according DIN EN ISO 20795-1).

Product shades: HinriDon is available in the following shades: colorless / clear, pink, pink opaque.

Storage and safety measures: The material must not be used after expiration date. Storage temperature must not exceed 25°C (77°F) Avoid direct sunlight. Keep containers closed after use. Keep away from heat / sparks / open flames / hot surfaces. No smoking. Dispose of contents/container in accordance with local recycling regulations. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid breathing dust / fume / gas / mist / vapours / spray. Wear protective gloves/protective clothing/eye protection/face protection. If on skin: Wash with plenty of soap and water

Contraindications: Patient contact with unpolymerized material or any of its components is contraindicated. If the patient is allergic to one or more ingredients of the product, the product must not be used.

Classification: Denture Base resin, short-term hot-curing acrylic based on methyl methacrylate in the form of powder and liquid. Classification according to ISO 20795-1 Type 1 Class 1 and according to medical device directive 93/42/EEC annex IX class IIa for removable dentures.

Distributor: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
Manufacturer: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH packaging carries the **CE** 0124 mark.

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El vendedor: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
El confeccionador: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH rotula el envase comercial con **CE** 0124.
Commerciale: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
Fabbicante: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar