# HinriPress®

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# HinriPress®

## Instructions for use

HinriPress® is the cold-curing all-rounder among all denture base acrylics and especially designed for dental technicians who like to cover all indications with only one product.

**Indications:**

- Total upper and lower dentures.
- Completion of model cast dentures
- Indirect relinings
- Partial upper and lower dentures
- Dilatations and repairs

**Technical data:**

Mixing ratio Powder/Liquid:	10 : 7
Thickening phase approx.	30 sec.*
Pouring phase approx.	3 min.*
Plastic phase approx.	3 min.*
Pressure polymerization at 45 - 50 °C approx.	7 min.*
Flashpoint:	10°C
Max. residual monomer content:	2,4 %

\* The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures extend, lower temperatures shorten these processing times.

**Product features:** cadmium free, very good flow ability, very high fitting accuracy, color stable due to a catalyst system without tertiary amines, mucosal friendly due to a low amount of residual monomer, fast and safe curing, certified biocompatibility.

Whether it is the completion of a total or partial denture, a simple repair, an indirect extraoral relining or a dilatation, HinriPress® always shows reliable results while saving precious time.

HinriPress® does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded), supplied with retentions and wetted with HinriPress®-Bonding.

**Preparation:** The arrangement of the prosthetic teeth in wax need to be fixed in a cuvette (HinriPress® – Cuvette G, HinriPress® – Cuvette S) or with a key made of plaster or silicone. Model and teeth are boiled out thoroughly to remove any wax residue. The plaster model needs to be well watered to assure that no air is pressed from the plaster into the acrylic while curing inside the pressure polymerization unit.

**Isolation:** Apply alginate based isolation on the model and keys of plaster. HinriPress K always shows reliable results. Important: The isolation needs to be well dried to make sure there are no white marks or spots on the finished acrylic.

**Processing:** Place required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as needed to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a thickening phase of approximately 30 sec. the material is pourable for another 2 to 3 minutes. Within this phase the material should be poured into the form. As soon as the plastic phase begins the material stops flowing and can be modelled for approx. 3 minutes.

**Polymerization:** After a maximum of 8 minutes the polymerization should be carried out in a pressure polymerization unit at a pressure of 2 to 3 bar and a water temperature of 45 to 50°C for 6 to 7 minutes.

**After completion:** Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste.

**Product shades:** HinriPress® is available in the following shades: colorless/clear, pink, pink opaque.

**Storage and safety measures:** Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapours. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

**Secondary effects:** Such as allergies to contents of the acrylic may very rarely occur.

**Reciprocal actions:** unknown

**Reactions:** In case of known incompatibilities and allergies to contents of the acrylic.

**Classification:** Denture base resin, Autopolymer based on methyl methacrylate consisting of powder and liquid. Classification according to ISO 20795-1 type 2 class 1 and according to medical device directive 93/42/EEC annex IX class IIa for removable and for fixed dentures.

Distributor: ERNST HINRICHS Dental GmbH, Borsigstraße 1, DE 38644 Goslar  
 Manufacturer: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar  
 SILADENT Dr. Böhme & Schöps GmbH packaging carries the 0124 mark.

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SILADENT Dr. Böhme & Schöps GmbH rotula el envase commercial con ☞ 0124.
Commerciale: ERNST HINRICHS Dental GmbH, Borsigstraße 1, DE 38644 Goslar
Fabbriente: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar