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Laser welding rods

Instructions for use

Area of application: Filler material for laser welding of CoCr dental alloys, which are suitable for laser welding.

Composition (in % by mass): Co: 65.30 Cr: 28.0 Mo: 5.5 **Other constituents:** Mn, Si

Contraindication: Do not use in cases of known allergy to any of the ingredients. Not to be used in combination with the customary flame or furnace soldering methods.

Precautions: When using a laser welder, the safety precautions provided by the manufacturer must be observed. To prevent any unnecessary material blending, the laser welding wires to be used should be of identical alloys.

Choice of material: Laser welding wires are used as laser welding filler material.

General recommendation for laser welding wires: Laser welding wires should not be used as solder with the usual flame or furnace soldering method.

Processing recommendations:

- Observe the manufacturer's recommendations for the device in question.
- Set the laser welder to the specific value required for the type of alloy used.
- Test the basic settings first on a piece of casting sprue whose surface has been pretreated with a silicone point and/or has been sandblast. Adjust the laser parameters as necessary.
- Caution: Reflective surfaces cannot be laser welded. The metallic area MUST be matt!
- Before laser welding, cast metal products should undergo heat treatment to prevent any material stresses
- Use of a protective gas prevents oxidation of the welding environment.
- Preferably use X-shaped laser welds.
- Build up the laser shots crosswise.
- We recommend that the object undergoes another heat. Treatment after laser welding to avoid material stress.

Safety hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask!

Further questions: Contact our technical services (Tel.: +49 (0) 53 21 – 5 06 24) or our sales representatives.

Distributor: ERNST HINRICHS Dental GmbH, Borsigstr. 1, DE 38644 Goslar
Manufacturer: EUKAMED OHG, Cranachstr. 44, 42549 Velbert

EUKAMED OHG packaging carries the mark

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El mercader: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Mercante: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar