Titan BioStar °5

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Titan BioStar °5

Instructions for use

Area of application: Non precious dental disc based on titanium, grade 5 for crowns and bridges technique and for metal-ceramic dentures accord. DIN EN ISO 22674, type 4 alloys Titan BioStar °5 is suitable for single crowns in anterior and posterior areas, frames in anterior and posterior areas with up to a 3 units span with adequate dimension.

Product description: Titan BioStar °5 is a titanium-based alloy for the dental application used for single crowns in front and lateral region, bridges in front and lateral region with up to 3 pontics with appropriate cross-section.

Composition (in % by mass)

Ti: ca. 90.0 Al: 5.5 – 6.75 V: 3.5 – 4.5

Other constituents:

Fe, O₂, C

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Technical data (guidelines)

Density g/cm ³	4.5
Elongation at rupture, minimum (%)	10
Tensile strength minimum (MPa):	895
Proof stress 0,2 % (MPa):	828
Modulus of elasticity (Gpa):	125
Melting point (°C):	1.700
WAK (25-500 °C):	10.2 x 10 ⁻⁶ K ⁻¹
Vickers hardness HV 10:	350
Casting temperature (°C)	1.660-1.700
Max. burning temperature (°C):	880

Grinding: Grind frames or individual abutments with tungsten carbide burs for titanium. Use also tungsten carbide burs or separating discs for cutting frames off the milling blank. All grinding can be done using cross-cut tungsten carbide burs or diamonds. It is important to cut always in the same direction, so as to avoid metal overlaps that might result in bubbles and porosities when firing the ceramics.

Attention: For top quality and efficiency, please follow all instructions closely. Always use protection gear (goggles, mask, suction...). Titan BioStar °5 dust and chips are highly inflammable! Therefore please follow safety precautions (no smoking when grinding, avoid open flames, etc.). TITAN BioStar °5 is only to be manipulated by trained staff.

Veneering: Use only titanium ceramics. Always follow ceramics manufacturer's instructions. Firing temperature should never exceed 880 °C, since a higher temperature would bring about molecular and structural alterations of the material. After veneering, rubber down all visible metal areas and polish using brushes and polishing pastes.

Soldering / Welding: To avoid oxide formation during the soldering / welding of titanium only use argon. To protect a corrosion resistance only use additional materials with the same ingredients.

Safety hint: Metal dust is hazardous to health. During the process use personal protective clothing (protecting glasses, dust mask, exhaust). For finishing and sandblasting use a suitable extraction system and/or face mask.

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

Storage conditions: No special storage conditions are required.

Further questions: Please contact ERNST HINRICHS technical services (Tel.: +49 (0) 53 21 – 5 06 24 / 25) or our sales representatives.

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Commerciale: ERNST HINRICHS GmbH, Borsigstr. 1, DE 38644 Goslar
Fabbricante: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH segnare il confezionamento a **CE 0124**
ERNST HINRICHS GmbH