HinriPress® Vario

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HinriPress® Vario

Instructions for use

HinriPress® Vario is a cold-curing all-rounder among all denture base acrylics. The expanded processing time of this material allows the dental technicians a most stress-free handling.

Technical data:

Mixing ratio Powder/Liquid:	10 : 7
Thickening phase approx.	30 sek.*
Pouring phase approx.	3 min.*
Plastic phase approx.	20 min.*
Pressure polymerization at 45 - 50 °C approx.	15 min.*
Flashpoint:	10 °C
Max. residual monomer content:	3,8 %

* The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures extend, lower temperatures shorten these processing times.

Indications:

- Completion of model cast dentures
- partial or total relinings (indirect method)
- Additions and repairs

Product features: Cadmium free, very long processing phase, very good flow ability, very high fitting accuracy, color stable due to a catalyst system without tertiary amines, mucosal friendly due to a low amount of residual monomer, fast and safe curing, certified biocompatibility. Whether it is the completion of a total or partial denture, a simple repair, an indirect extraoral relining or a dilatation, HinriPress® Vario always shows reliable results. **HinriPress® Vario does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded), supplied with retentions and wetted with HinriPress®-Bonding.**

Preparation: The arrangement of the prosthetic teeth in wax need to be fixed with a key made of gypsum (e. G. Hinridur) or silicone. Model and teeth are boiled out thoroughly to remove any wax residue. The gypsum model needs to be well watered to assure that no air is pressed from the gypsum into the acrylic while curing inside the pressure polymerization unit.

Isolation: Apply alginate based isolation on the model and keys of plaster. HinriSep K always shows reliable results. **Important:** The isolation needs to be well dried to make sure there are no white marks or spots on the finished acrylic.

Processing: Place required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as needed to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a thickening phase of approximately 30 sec. the material is pourable for another 3 minutes. Within this phase the material should be poured into the form. As soon as the plastic phase begins the material stops flowing and can be modelled for approx. 20 minutes.

Polymerization: After a maximum of 20 minutes the polymerization should be carried out in a pressure polymerization unit at a pressure of 2 to 3 bar and a water temperature of 45 to 50°C for approx. 15 minutes.

After completion: Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste (3.8 % according DIN EN ISO 20795-1).

Product shades: HinriPress® Vario is available in the following shades: colorless / clear, pink, pink opaque.

Storage and safety measures: The material must not be used after expiration date. Storage temperature must not exceed 25°C (77°F) Avoid direct sunlight. Keep containers closed after use. Keep away from heat / sparks / open flames / hot surfaces. No smoking. Dispose of contents/container in accordance with local recycling regulations. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid breathing dust / fume / gas / mist / vapours / spray. Wear protective gloves/protective clothing/eye protection/face protection. If on skin: Wash with plenty of soap and water

Contraindications: Patient contact with unpolymerized material or any of its components is contraindicated. If the patient is allergic to one or more ingredients of the product, the product must not be used.

Reactions: In case of known incompatibilities and allergies to contents of the acrylic.

Classification: Denture base resin, Autopolymer based on methyl methacrylate consisting of powder and liquid. Classification according to ISO 20795-1 type 2 class 1 and according to medical device directive 93/42/EEC annex IX class IIa for removable and for fixed dentures.

Distributor: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
Manufacturer: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH packaging carries the **CE 0124** mark.

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El vendedor: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
El confeccionador: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH rotula el envase comercial con **CE** 0124.
Commerciale: ERNST HINRICHS GmbH, Borsigstraße 1, DE 38644 Goslar
Fabbricante: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE-38644 Goslar