




according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N


SECTION 1: Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
 Tradename: Hinri-Alloy N
 Product description: Nickel alloy
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
 Relevant identified uses: Manufacturing of dental prosthesis in dental laboratories
- 1.3 Details of the supplier of the safety data sheet
 Manufacturer/Supplier: ERNST HINRICHS Dental GmbH
 Street / mailbox: Borsigstr. 1
 Country code. / postal code / city: D - 38644 Goslar
 Phone: 0 53 21 / 5 06 24
 Fax: 0 53 21 / 5 08 81
 E-mail / Website: info@hinrichs-dental.de / www.hinrichs-dental.de
 Further information obtainable from: ERNST HINRICHS Dental GmbH
- 1.4 Emergency telephone number
 ERNST HINRICHS Dental GmbH: +49 (0) 53 21 / 5 06 24 - 25 (Mon-Fri. 8 a.m. - 4 p.m.)

SECTION 2: Hazards identification

- 2.1 Classification of the substance or mixture
 Classification according to Regulation (EC) No 1272/2008:
 Skin Sens. 1 H317 May cause an allergic skin reaction.
 Resp. Sens. 1 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
 Aquatic Chronic 4 H413 May cause long lasting harmful effects to aquatic life
 Adverse human health and environmental effects: No additional information available.
- 2.2 Label elements
 Labelling according to Regulation (EC) No 1272/2008:
 Labelling is not applicable according to Annex I; 1.3.4. (Metals in massive form, alloys, mixtures containing polymers, mixtures containing elastomers).
 In the form in which it is marketed, the product causes no danger to health for humans through inhalation, swallowing or contact with the skin.
 Processing by downstream users causing differences in the product form for which Regulation (EC) No. 1272/2008 Annex I; 1.3.4 does not apply are obliged to label.
 Hazard pictograms:
- 

GHS07



GHS08
- Signal word: Danger
 Hazard-determining components of labelling: Nickel
- Hazard statements:
 H317 May cause an allergic skin reaction.
 H351 Suspected of causing cancer.
 H372 Causes damage to organs through prolonged or repeated exposure.
- Precautionary statements:
 P201 Obtain special instructions before use.
 P308 + P313 IF exposed: Call a POISON CENTER or doctor/physician.
 P314 Get Medical advice/attention if you feel unwell.
 P363 Wash contaminated clothing before reuse.
 P405 Store locked up.
- 2.3 Other hazards
 Results of PBT and vPvB assessment:
 PBT: Not applicable.
 vPvB: Not applicable.



according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N

SECTION 3: Composition/information on ingredients

3.1 Substances Not applicable

3.2 Mixtures
 Description: alloy
 Dangerous components:

CASAS-No.7440-02-0 EC-No. 231-111-4	Nickel	63 %
	Carc. 2; H351 Skin Sens. 1; H317 STOT RE 1; H372	
CAS: 7440-47-3 EINECS: 231-157-5	Chromium	25 %
	Not classified	
	Not classified	

Additional information: For the wording of the listed risk phrases refer to section 16.

SECTION 4: First aid measures

- 4.1 Description of first aid measures
 General information: Instantly remove any clothing soiled by the product.
 Take affected persons into the open air.
 After inhalation: In case of unconsciousness bring patient into stable side position for transport. Supply fresh air.
 After skin contact: Instantly wash with water and soap and rinse thoroughly. If skin irritation continues, consult a doctor.
 After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult doctor.
 After swallowing: Rinse mouth. Drink plenty of water as a precaution. Get medical advice/attention.
- 4.2 Most important symptoms and effects, both acute and delayed: No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed: No further relevant information available.

SECTION 5: Firefighting measures

- 5.1 Extinguishing media
 Suitable extinguishing agents: Metal fire powders, Sand
 For safety reasons unsuitable extinguishing agents: Water, Foam, CO₂
- 5.2 Special hazards arising from the substance or mixture: No additional information available.
- 5.3 Advice for firefighters
 Protective equipment: Do not inhale explosion gases or combustion gases.
 Wear self-contained breathing apparatus.
 Additional information: Dispose of fire debris and contaminated firefighting water in accordance with official regulations.
 Collect contaminated firefighting water separately. It must not enter drains.

SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures: Wear protective equipment. Keep unprotected persons away.
 Ensure adequate ventilation
 Use breathing protection against the effects of dust.
- 6.2 Environmental precautions: Prevent material from reaching sewage system, holes and cellars.
 Inform respective authorities in case product reaches water or sewage system.
 Do not allow to enter drainage system, surface or ground water.



according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N

- 6.3 Methods and material for containment and cleaning up: Take up mechanically and collect in suitable container for disposal. Send for recovery or disposal in suitable containers.
- 6.4 Reference to other sections: See Section 7 for information on safe handling
 See Section 8 for information on personal protection equipment.
 See Section 13 for information on disposal.

SECTION 7: Handling and storage

- 7.1 Precautions for safe handling: Store in cool, dry place in tightly closed containers. Ensure good ventilation/exhaustion at the workplace. Take note of emission threshold. Keep away from children. Avoid release to the environment.
 Information about protection against explosions and fires: Dust can form an explosive mixture with air. Avoid formation of dust.
- 7.2 Conditions for safe storage, including any incompatibilities: Storage: No further relevant information available.
 Storage class: 10-13
- 7.3 Specific end use(s): No further relevant information available.

SECTION 8: Exposure controls/personal protection

- 8.1 Control parameters:

Components with critical values that require monitoring at the workplace:	
7440-02-0 Nickel	
WEL	Long-term value: 0.1 mg/m ³ Skin resorption/sensibilisation: Sk
7440-47-3 Chromium	
WEL	Long-term value: 0,5 mg/m ³
- 8.2 Exposure controls:
 - Additional information: Based on information valid at the time of writing.
 - Additional information about design of technical systems: Please take care on national and local requirements.
 - General protective and hygienic measures: Keep away from food, drink and animal feedingstuffs. Instantly remove any soiled and impregnated garments. Wash hands during breaks and at the end of the work. Avoid contact with the eyes and skin.
 - Personal protective equipment:
 - Protection of hands: Wear suitable gloves. The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
 - Material of gloves: Latex, Nitrile rubber (> 0.1 mm)
 - Penetration time of glove material: The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
 - Eye protection: Tightly sealed safety glasses (EN 166).
 - Breathing equipment: In case of dust formation or inadequate ventilation wear respiratory protection. Breathing apparatus with filter P. Not necessary if room is well-ventilated.
 - Body protection: Protective work clothing.

SECTION 9: Physical and chemical properties

- 9.1 Information on basic physical and chemical properties:
 - General Information
 - Appearance: Solid



according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N

Physical state:	Solid
Colour:	metallic
Odour:	No data available
Odour threshold:	No data available
pH:	No data available
Melting point/Freezing point	1380 °C
Boiling point/Boiling range:	No data available
Flash point:	No data available
Flammability (solid, gas):	No data available.
Auto-ignition temperature:	Product is not self-igniting.
Decomposition temperature:	No data available
Oxidising properties	No data available
Explosive properties:	Product is not explosive. However, during processing formation of explosive dusts is possible.
Explosive Limits:	
Lower:	No data available
Upper:	No data available
Vapour pressure at 20 °C:	No data available
Density at 20 °C:	8,3 g/cm ³
Relative density:	No data available
Vapour density:	No data available.
Evaporation rate:	No data available
Solubility:	Insoluble.
Partition coefficient (n-octanol/water):	No data available
Viscosity:	
dynamic:	No data available
kinematic:	No data available
9.2 Other information	No further relevant information available.

SECTION 10: Stability and reactivity

10.1 Reactivity	No dangerous reactions known.
10.2 Chemical stability	
Thermal decomposition / conditions to be avoided:	No decomposition if used according to specifications.
10.3 Possibility of hazardous reactions	Possible hydrogen formation upon contact with acids.
10.4 Conditions to avoid	None, if handled according to intended use.
10.5 Incompatible materials	Acids.
10.6 Hazardous decomposition products	No dangerous decomposition products known.

SECTION 11: Toxicological information

11.1 Information on toxicological effects	
Acute toxicity:	Not classified Based on available data, the classification criteria are not met.
Skin corrosion/irritation:	Not classified Based on available data, the classification criteria are not met.
Serious eye damage/irritation:	Not classified Based on available data, the classification criteria are not met.
Respiratory or skin sensitisation:	May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction.
Germ cell mutagenicity:	Not classified Based on available data, the classification criteria are not met.



according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N

Carcinogenicity:	Not classified Based on available data, the classification criteria are not met.
Reproductive toxicity:	Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (single exposure):	Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (repeated exposure):	Not classified Based on available data, the classification criteria are not met.
Aspiration hazard:	Not classified Based on available data, the classification criteria are not met.
Potential adverse human health effects and symptoms:	Inhalation of vapours in high concentration may cause metal fume fever and may cause damage of the central nervous system in case repeated and prolonged exposure. Repeated and prolonged exposure to high dust concentrations may lead to irritation of the respiratory tract. Inhalation of metal-containing dusts may cause acute poisoning, leading to nausea, vomiting and abdominal pain.

SECTION 12: Ecological information

12.1 Toxicity Aquatic toxicity:	May cause long lasting harmful effects to aquatic life.
12.2 Persistence and degradability:	No further relevant information available.
12.3 Bioaccumulative potential:	No further relevant information available.
12.4 Mobility in soil:	No further relevant information available.
12.5 Results of PBT and vPvB assessment PBT: vPvB:	Not applicable. Not applicable.
12.6 Other adverse effects	No further relevant information available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods Recommendation:	Disposal must be made according to official regulations. Do not dispose of with domestic waste. Do not allow product to reach sewage system.
---	--

European waste catalogue	
18 00 00	WASTES FROM HUMAN OR ANIMAL HEALTH CARE AND/OR RELATED RESEARCH (EXCEPT KITCHEN AND RESTAURANT WASTES NOT ARISING FROM IMMEDIATE HEALTH CARE)
18 01 00	wastes from natal care, diagnosis, treatment or prevention of disease in humans
18 01 06*	chemicals consisting of or containing dangerous substances

SECTION 14: Transport information

14.1 UN number ADR, IMDG, IATA:	Not regulated for transport
14.2 UN proper shipping name ADR, IMDG, IATA:	Not applicable
14.3 Transport hazard class(es) ADR, IMDG, IATA Klasse:	Not applicable



according 1907/2006/EG, Article 31
 Date of issue: 03.04.2013
 Date of revision: 11.01.2016 / Version: 3.0

Hinri-Alloy N

Danger labels:	Not applicable
14.4 Packing group ADR, IMDG, IATA:	Not applicable
14.5 Environmental hazards	No
Marine pollutant:	No
Special marking:	No supplementary information available
14.6 Special precautions for user	Not applicable
Kemler Number:	Not applicable
EMS Number:	Not applicable
14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture	
National regulations:	
Other regulations, limitations and prohibitive regulations:	
15.2 Chemical safety assessment	A Chemical Safety Assessment of the product has not been carried out.

SECTION 16: Other information

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. The data given here only applies when product used for proper application(s). The product is not sold as suitable for other applications - usage in such may cause risks not mentioned in this sheet. Do not use for other application(s) without seeking advice from manufacturer.

Adjustments compared to earlier versions
 Section 1-16 New edition, complete revision

Relevant phrases

H317:	May cause an allergic skin reaction
H334:	May cause allergy or asthma symptoms or breathing difficulties if inhaled
H413:	May cause long lasting harmful effects to aquatic life
H351	Suspected of causing cancer.
H372	Causes damage to organs through prolonged or repeated exposure.

Abbreviations and acronyms:

ADR:	Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG:	International Maritime Code for Dangerous Goods
IATA:	International Air Transport Association
GHS:	Globally Harmonised System of Classification and Labelling of Chemicals
EINECS:	European Inventory of Existing Commercial Chemical Substances
CAS:	Chemical Abstracts Service (division of the American Chemical Society)
Aquatic Chronic 4:	Hazardous to the aquatic environment — Chronic Hazard, Category 4
Resp. Sens. 1:	Sensitisation — Respiratory, category 1
Skin Sens. 1:	Sensitisation — Skin, category 1