



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Keralloy BioStar

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
 Tradename: Keralloy BioStar
 Product description: Cobalt-chromium alloy
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
 Relevant identified uses: Manufacturing of dental prosthesis in dental laboratories
- 1.3 Details of the supplier of the safety data sheet
 Manufacturer/Supplier: ERNST HINRICHS Dental GmbH
 Street / mailbox: Borsigstr. 1
 Country code. / postal code / city: D - 38644 Goslar
 Phone: 0 53 21 / 5 06 24
 Fax: 0 53 21 / 5 08 81
 E-mail / Website: info@hinrichs-dental.de / www.hinrichs-dental.de
 Further information obtainable from: ERNST HINRICHS Dental GmbH
- 1.4 Emergency telephone number
 ERNST HINRICHS Dental GmbH: +49 (0) 53 21 / 5 06 24 - 25 (Mon-Fri. 8 a.m. – 4 p.m.)

SECTION 2: Hazards identification

- 2.1 Classification of the substance or mixture
 Classification according to Regulation (EC) No 1272/2008:
 Skin Sens. 1 H317 May cause an allergic skin reaction.
 Resp. Sens. 1 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
 Aquatic Chronic 4 H413 May cause long lasting harmful effects to aquatic life
 Adverse human health and environmental effects: No additional information available.
- 2.2 Label elements
 Labelling according to Regulation (EC) No 1272/2008:
 Labelling is not applicable according to Annex I; 1.3.4. (Metals in massive form, alloys, mixtures containing polymers, mixtures containing elastomers).
 In the form in which it is marketed, the product causes no danger to health for humans through inhalation, swallowing or contact with the skin.
 Processing by downstream users causing differences in the product form for which Regulation (EC) No. 1272/2008 Annex I; 1.3.4 does not apply are obliged to label.
 Hazard pictograms:
- GHS08
 Danger
 Cobalt
- Signal word: Danger
 Hazard-determining components of labelling: Cobalt
- Hazard statements:
 H317 May cause an allergic skin reaction.
 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
 H413 May cause long lasting harmful effects to aquatic life.
- Precautionary statements:
 P273: Avoid release to the environment.
 P284: In case of inadequate ventilation wear respiratory protection.
 P302+P352: IF ON SKIN: Wash with plenty of water and soap.
 P304+P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing.
 P342+P311: If experiencing respiratory symptoms: Call a doctor.
 P362+P364: Take off contaminated clothing and wash it before reuse.
- 2.3 Other hazards
 Results of PBT and vPvB assessment:
 PBT: Not applicable.



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Page 2 of 7
 Printing date: 16.02.2018

Keralloy BioStar

vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

3.1 Substances Not applicable

3.2 Mixtures
 Description: Alloy
 Dangerous components:

CAS: 7440-48-4 EINECS: 231-158-0 Index: 027-001-00-9	Cobalt	66%
	Skin Sens. 1, H317; Resp. Sens. 1, H334; Aquatic Chronic 4, H413	
CAS: 7440-47-3 EINECS: 231-157-5	Chromium	27%
	Not classified	
CAS: 7439-98-7 EINECS: 231-107-2	Molybdenum	6%
	Not classified	
CAS: 7440-21-3 EINECS: 215-609-9	Silicon	< 1%
	Not classified	
CAS: 7439-96-5 EINECS: 231-105-1	Manganese	< 1%
	Not classified	

Additional information: For the wording of the listed risk phrases refer to section 16.

SECTION 4: First aid measures

- 4.1 Description of first aid measures
 General information: Instantly remove any clothing soiled by the product. Take affected persons into the open air.
 After inhalation: In case of unconsciousness bring patient into stable side position for transport. Supply fresh air.
 After skin contact: Instantly wash with water and soap and rinse thoroughly. If skin irritation continues, consult a doctor.
 After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult doctor.
 After swallowing: Rinse mouth. Drink plenty of water as a precaution. Get medical advice/attention.
- 4.2 Most important symptoms and effects, both acute and delayed: May cause an allergic skin reaction. May cause allergy or asthma symptoms or breathing difficulties if inhaled.
- 4.3 Indication of any immediate medical attention and special treatment needed: Treat symptomatically.

SECTION 5: Firefighting measures

- 5.1 Extinguishing media
 Suitable extinguishing agents: Metal fire powders, Sand
 For safety reasons unsuitable extinguishing agents: Water, Foam, CO₂
- 5.2 Special hazards arising from the substance or mixture: No additional information available.
- 5.3 Advice for firefighters
 Protective equipment: Do not inhale explosion gases or combustion gases. Wear self-contained breathing apparatus.
 Additional information: Dispose of fire debris and contaminated firefighting water in accordance with official regulations. Collect contaminated firefighting water separately. It must not enter drains.



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Keralloy BioStar

SECTION 6: Accidental release measures

- | | | |
|-----|--|---|
| 6.1 | Personal precautions, protective equipment and emergency procedures: | Wear protective equipment. Keep unprotected persons away. Ensure adequate ventilation. Use breathing protection against the effects of dust. |
| 6.2 | Environmental precautions: | Do not allow to enter drainage system, surface or ground water. Inform respective authorities in case product reaches water or sewage system. |
| 6.3 | Methods and material for containment and cleaning up: | Take up mechanically and collect in suitable container for disposal. Send for recovery or disposal in suitable containers. |
| 6.4 | Reference to other sections: | See section 7 for information on safe handling.
See section 8 for information on personal protection equipment.
See section 13 for information on disposal. |

SECTION 7: Handling and storage

- | | | |
|-----|--|---|
| 7.1 | Precautions for safe handling

Information about protection against explosions and fires:
Hygiene measures: | Ensure good ventilation/exhaustion at the workplace. Do not breathe dust.

Dust can form an explosive mixture with air. Avoid formation of dust.

When handling the product, ensure good industrial hygiene and adequate safety precautions. Do not eat, drink or smoke at work. Before eating, drinking, smoking and leaving the workplace, wash hands and other bare areas with mild soap and water. Do not wear contaminated clothing outside the workplace. |
| 7.2 | Conditions for safe storage, including any incompatibilities
Storage conditions:
Information on mixed storage:

Storage class: | Store in cool, dry place in tightly closed containers.
Keep away from food, drink and animal feedingstuffs. Do not store together with explosive substances.
10-13 |
| 7.3 | Specific end use(s) | No further relevant information available. |

SECTION 8: Exposure controls/personal protection

- | | | |
|---|---|--|
| 8.1 | Control parameters | |
| Components with critical values that require monitoring at the workplace: | | |
| 7440-48-4 Cobalt | | |
| WEL | Long-term value: 0.1 mg/m ³ | |
| 7440-47-3 Chromium | | |
| WEL | Long-term value: 0.5 mg/m ³ | |
| 7439-98-7 Molybdenum | | |
| WEL | Long-term value: 10 mg/m ³
Short-term value: 20 mg/m ³ | |
| 7440-21-3 Silicon | | |
| WEL | Long-term value: 10 mg/m ³ inhalable dust, 4 mg/m ³ respirable dust | |
| 7439-96-5 Manganese | | |
| WEL | Long-term value: 0.5 mg/m ³ | |
| 8.2 | Exposure controls:
Appropriate engineering controls: | Ensure adequate ventilation so that concentrations do not exceed the applicable standard values. |



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Keralloy BioStar

Personal protective equipment:	
Protection of hands:	Wear suitable gloves. The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Material of gloves:	Latex, Nitrile rubber (> 0.1 mm).
Penetration time of glove material:	The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
Eye protection:	Tightly sealed safety glasses (EN 166).
Breathing equipment:	In case of dust formation or inadequate ventilation wear respiratory protection. Breathing apparatus with filter P. Not necessary if room is well-ventilated.
Body protection:	Protective work clothing.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

General Information

Appearance:	Solid
Physical state:	Solid
Colour:	Metallic
Odour:	No data available
Odour threshold:	No data available
pH:	No data available
Melting point/Freezing point	1436 °C
Boiling point/Boiling range:	No data available
Flash point:	Not applicable
Flammability (solid, gas):	No data available.
Auto-ignition temperature:	Product is not self-igniting.
Decomposition temperature:	No data available
Oxidising properties	No data available
Explosive properties:	Product is not explosive. However, during processing formation of explosive dusts is possible.
Explosive Limits:	
Lower:	No data available
Upper:	No data available
Vapour pressure at 20 °C:	No data available
Density at 20 °C:	8.5 g/cm ³
Relative density:	No data available
Vapour density:	No data available
Evaporation rate:	No data available
Solubility:	Insoluble
Partition coefficient (n-octanol/water):	No data available
Viscosity:	
dynamic:	No data available
kinematic:	No data available

9.2 Other information

No further relevant information available.

SECTION 10: Stability and reactivity

10.1 Reactivity	Under normal conditions of use no dangerous reactions known.
10.2 Chemical stability	
Thermal decomposition / conditions to be avoided:	Stable under recommended storage and application conditions according to section 7. No decomposition if used according to specifications.



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Keralloy BioStar

10.3 Possibility of hazardous reactions	Possible hydrogen formation upon contact with acids.
10.4 Conditions to avoid	None, if handled according to intended use.
10.5 Incompatible materials	Acids.
10.6 Hazardous decomposition products	No dangerous decomposition products known.

SECTION 11: Toxicological information

11.1 Information on toxicological effects	
Acute toxicity:	Not classified Based on available data, the classification criteria are not met.
Skin corrosion/irritation:	Not classified Based on available data, the classification criteria are not met.
Serious eye damage/irritation:	Not classified Based on available data, the classification criteria are not met.
Respiratory or skin sensitisation:	May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction.
Germ cell mutagenicity:	Not classified Based on available data, the classification criteria are not met.
Carcinogenicity:	Not classified Based on available data, the classification criteria are not met.
Reproductive toxicity:	Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (single exposure):	Not classified Based on available data, the classification criteria are not met.
Specific target organ toxicity (repeated exposure):	Not classified Based on available data, the classification criteria are not met.
Aspiration hazard:	Not classified Based on available data, the classification criteria are not met.
Potential adverse human health effects and symptoms:	Inhalation of vapours in high concentration may cause metal fume fever and may cause damage of the central nervous system in case repeated and prolonged exposure. Repeated and prolonged exposure to high dust concentrations may lead to irritation of the respiratory tract. Inhalation of metal-containing dusts may cause acute poisoning, leading to nausea, vomiting and abdominal pain.

SECTION 12: Ecological information

12.1 Toxicity	
Aquatic toxicity:	May cause long lasting harmful effects to aquatic life.
12.2 Persistence and degradability:	No further relevant information available.
12.3 Bioaccumulative potential:	No further relevant information available.
12.4 Mobility in soil:	No further relevant information available.
12.5 Results of PBT and vPvB assessment	
PBT:	Not applicable.
vPvB:	Not applicable.
12.6 Other adverse effects	No further relevant information available.



according to Regulation (EU) 2015/830
 Date of issue: 03.04.2013
 Date of revision: 15.02.2018 / Version: 4.0

Page 6 of 7
 Printing date: 16.02.2018

Keralloy BioStar

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Recommendation: Disposal must be made according to official regulations. Do not dispose of with domestic waste. Do not allow product to reach sewage system.

European waste catalogue	
18 00 00	WASTES FROM HUMAN OR ANIMAL HEALTH CARE AND/OR RELATED RESEARCH (EXCEPT KITCHEN AND RESTAURANT WASTES NOT ARISING FROM IMMEDIATE HEALTH CARE)
18 01 00	wastes from natal care, diagnosis, treatment or prevention of disease in humans
18 01 06*	chemicals consisting of or containing dangerous substances

Waste code: The waste code number according to the Ordinance on the European Waste Catalogue (AVV) depends on the waste producer and can therefore vary for any given product. The waste code number is therefore to be gleaned separately from each waste producer.

SECTION 14: Transport information

- 14.1 UN number
 ADR, IMDG, IATA: Not regulated for transport
- 14.2 UN proper shipping name
 ADR, IMDG, IATA: Not applicable
- 14.3 Transport hazard class(es)
 ADR, IMDG, IATA
 Class: Not applicable
 Danger labels: Not applicable
- 14.4 Packing group
 ADR, IMDG, IATA: Not applicable
- 14.5 Environmental hazards
 Marine pollutant: No
 Special marking: No supplementary information available
- 14.6 Special precautions for user
 Kemler Number: Not applicable
 EMS Number: Not applicable
- 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code: Not applicable

SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
 National regulations:
 Other regulations, limitations and prohibitive regulations:
- 15.2 Chemical safety assessment: A Chemical Safety Assessment of the product has not been carried out.

SECTION 16: Other information

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. The data given here only applies when product used for proper application(s). The product is not sold as suitable for other applications - usage in such may cause risks not mentioned in this sheet. Do not use for other application(s) without seeking advice from manufacturer.
 Adjustments compared to earlier versions
 Section 1-16 New edition, complete revision



according to Regulation (EU) 2015/830

Date of issue: 03.04.2013

Date of revision: 15.02.2018 / Version: 4.0

Keralloy BioStar

Relevant phrases

H317: May cause an allergic skin reaction

H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled

H413: May cause long lasting harmful effects to aquatic life

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route
(European Agreement concerning the International Carriage of Dangerous Goods
by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

Aquatic Chronic 4: Hazardous to the aquatic environment — Chronic Hazard, Category 4

Resp. Sens. 1: Sensitisation — Respiratory, category 1

Skin Sens. 1: Sensitisation — Skin, category 1