



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 1 of 6

Printing date: 19.05.2016

**Wax BioStar**

**1. Identification of the substance / Preparation and Company:**

|     |   |   |
|-----|---|---|
| 1.1 | Product identifier  |   |
|     | Tradename:  | Wax BioStar   |
| 1.2 | Relevant identified uses of the substance or mixture and uses advised against |   |
|     | Relevant identified uses:   | to be used in dental labs exclusively                   |
|     | Uses advised against:   | should not be used by untrained personnel               |
| 1.3 | Details of the supplier of the safety data sheet                              |   |
|     | Manufacturer/Supplier:  | ERNST HINRICHS Dental GmbH                              |
|     | Street / mailbox:   | Borsigstr. 1  |
|     | Country code. / postal code / city:   | D - 38644 Goslar  |
|     | Phone:  | 0 53 21 / 5 06 24                                       |
|     | Fax:  | 0 53 21 / 5 08 81                                       |
|     | E-mail / Website:   | info@hinrichs-dental.de / www.hinrichs-dental.de        |
|     | Further information obtainable from:  | ERNST HINRICHS Dental GmbH                              |
| 1.4 | Emergency telephone number  |   |
|     | ERNST HINRICHS Dental GmbH:   | +49 (0) 53 21 / 5 06 24 - 25 (Mon-Fri. 8 a.m. – 4 p.m.) |

**2. Hazards identification:**

|     |  |  |
|-----|--|--|
| 2.1 | Classification of the substance or mixture:                    | Classification according to regulation (EC) no. 1272/2008 (CLP)  |
|     |  | This substance does not meet the criteria for classification according to regulation (EC) no. 1272/2008 (CLP). |
| 2.2 | Label Elements   |  |
|     | Labelling according to regulation (EC) no. 1272/2008 (CLP/GHS) |  |
|     | Hazard pictograms:   | NONE   |
|     | Signal word:   | NONE   |
|     | Hazardous component(s) for labelling:                          |  |
|     | Hazard statements:   | NONE   |
|     | Precautionary statements:                                      |  |
|     | P210:  | Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.                 |
|     | P262:  | Do not get in eyes, on skin, or on clothing.   |
|     | P370+  | In case of fire:   |
|     | P378:  | Use foam, carbon dioxide, dry powder to extinguish.  |
|     | Supplemental Hazard information (EU):                          | NONE   |
| 2.3 | Other Hazards:   | Results of PBT and vPvP assessment:  |
|     |  | PBT not applicable   |
|     |  | vPvP not applicable  |

**3. Composition/information on ingredients:**

|     |                                       |                         |
|-----|---------------------------------------|-------------------------|
| 3.1 | Substances                            |                         |
|     | Principal component of the substance: |                         |
|     | Substance name:                       | paraffin                |
|     | Index-no.:                            | void                    |
|     | EG-no.:                               | 232-315-6               |
|     | CAS-no.:                              | 8002-74-2               |
|     | Impurities, stabilizers:              | none                    |
|     | Substance name:                       |                         |
|     | Index-no.:                            |                         |
|     | EG-no.:                               |                         |
|     | CAS-no.:                              |                         |
| 3.2 | Mixtures                              | none                    |
|     | Description of the mixture:           |                         |
|     | EG-no.:                               | CAS-no. :               |
|     |                                       | Index-Nno.:             |
|     |                                       | REACH-registration no.: |



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 2 of 6

Printing date: 19.05.2016

**Wax BioStar**

Content ration: %  
 Classification according to regulation (EC) none  
 no. 1272/2008 (CLP):  
 (Full text of H- and EUH-phrases: see SECTION 16)

**4. First aid measures:**

- 4.1 Description of first aid measures
- Following inhalation: Pure wax resp. wax not containing any kind of oil is showing almost no vapour pressure when fluid, for that an inhalative absorption of any toxicologically concerned component can be excluded.
- Following skin contact: Contact with molten/fluid product may lead to burns. Cool down with cold water immediately, if necessary treat wound in a sterile manner, seek for medical advice. Any absorption of paraffin components by means of skin while the wax is solid cannot be expected.
- Following eye contact: Rinse thoroughly with water, if necessary seek medical advice (when wax is hot/fluid)
- Following ingestion: Rinse mouth – no further measures necessary. Waxes can be applied as components of pharmaceuticals or food for this reason the question of absorption in the digestive tract after accidental ingestion is irrelevant.
- 4.2 Most important symptoms and effects, both acute and delayed  
 Symptoms: none  
 Effects: none
- 4.3 Indication of any immediate medical attention and special treatment needed  
 Notes for the doctor: none  
 Special treatment: none

**5. Firefighting measures:**

- 5.1 Extinguishing media  
 Suitable: Foam, dry powder, carbon dioxide, water mist spray  
 Unsuitable: Water jet
- 5.2 Special hazards arising from the substance or mixture: In case of fire fumes are possible: carbon monoxide (CO).
- 5.3 Advice for fire-fighters: Respiratory protection independent of ambient environment – in case of fumes.
- 5.4 Additional indication:  
 Fire class: B: Liquid substances or substances becoming liquid.

**6. Accidental release measures:**

- 6.1 Personal precautions, protective equipment and emergency procedures  
 For non-emergency personnel  
 Protective equipment:  
 Emergency procedures:  
 Avoid contact with skin, eyes and clothes when substance is hot and fluid.  
 For emergency responders  
 Personal protective equipment: see section 5.3  
 Avoid contact with skin, eyes and clothes when substance is hot and fluid.
- 6.2 Environmental precautions: No special precautions need to be taken.
- 6.3 Methods and material for containment and cleaning up: Let material solidify – clean up mechanically.
- 6.4 Reference to other sections: More information on disposal see section 13.



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 3 of 6

Printing date: 19.05.2016

**Wax BioStar**

**7. Handling and storage:**

- 7.1 Precautions for safe handling
- Fire preventions: No special precautions necessary. Regular preventive fire protection.
- Aerosol and dust generation preventions: Avoid generation of dust – work under fume hood.
- Environmental precautions: No special precautions necessary.
- Advice on general occupational hygiene: Observe general hygiene rules.
- 7.2 Conditions for safe storage, including any incompatibilities
- Storage conditions: Keep containers tightly closed and store cool, dry and well ventilated. Keep away from sun rays. Best at 5 – 30 °C. Store in original bundles. Room temperature see above.
- Requirements for storage rooms and vessels: 11 (TRGS 510) Materials to avoid: n.a.
- Storage class: 11 (TRGS 510) Materials to avoid: n.a.
- 7.3 Specific end uses
- Sector specific guidelines: dental laboratories

**8. Exposure controls / personal protection:**

- 8.1 Control parameters
- Occupational exposure limits
- Substance: paraffin: CAS-no. : 8002-74-2
- Specification: Does not contain any substances - according to the common EU lists – that have to be monitored at work-places.
- Limit value type (country of origin): n.a.
- Peak limitation: n.a.
- Toxic to reproduction: n.a.
- Source (TRGS): n.a.
- 8.1.2 DNEL- and PNEC- values: none
- Substance: CAS-Nr. :
- Specification:
- Value:
- 8.1.3 Control-Banding (e.g. ILO, EMKG): none
- 8.2 Exposure controls:
- The type and extent of usage determine the choice of protective measures.
- Appropriate engineering controls: Provide good ventilation
- Personal protective equipment:
- Eye / face protection: Safety glasses – not necessary if used according to the intended purpose.
- Skin protection
- Hand protection
- in full contact: use appropriate safety gloves
- material: polychloroprene - CR (0,5 mm), nitrile rubber - NBR (0,35 mm)
- layer thickness (mm): see above
- penetration time (min.): 8 hours
- Splash contact: use appropriate safety gloves
- material: polychloroprene - CR (0,5 mm), nitrile rubber - NBR (0,35 mm)
- layer thickness (mm): see above
- penetration time (min.): 8 hours
- Body protection: Appropriate
- Respiratory protection: Not required under normal use.
- Thermal hazards: Not required under normal use.
- Environmental exposure controls: Do not allow to get into sewer system, surface or ground water.



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 4 of 6

Printing date: 19.05.2016

**Wax BioStar**

**9. Physical and chemical properties:**

9.1 Information on basic physical and chemical properties:

|                                       |   |
|---------------------------------------|---|
| Physical state:                       | Solid   |
| Aggregate state:                      | Ivory   |
| colour :                              | Characteristic  |
| Odour :                               | No information available                                  |
| Odour threshold :                     | n.a.  |
| pH-value:                             | < 50°C  |
| Freezing point :                      | n.a.  |
| Initial boiling point/boiling range : | > 130°C   |
| Flash point:                          | n.a.  |
| Evaporation rate :                    | The substance is flammable, but does not intensify fires. |
| Flammability:                         | none  |
| Upper/lower explosive limits          | n.a.  |
| Vapour pressure :                     | n.a.  |
| Vapour density:                       | n.a.  |
| Relative density:                     | ca. 0,93 g/cm <sup>3</sup>                                |
| Solubility:                           | Not soluble in water                                      |
| Partition coefficient:                | n.a.  |
| n-octanol/water:                      |   |
| Auto-ignition temperature:            | n.a.  |
| Decomposition temperature:            | n.a.  |
| Viscosity:                            | n.a.  |
| explosive properties:                 | none  |
| Oxidising properties:                 | none  |

9.2 Other information

Physical Explosives hazards: No further information available

**10. Stability and reactivity:**

|  |                                  |
|--|----------------------------------|
| 10.1 Reactivity:                         | No known unusual reactivity.     |
| 10.2 Chemical stability:                 | Stable under regular conditions. |
| 10.3 Possibility of hazardous reactions: | No known hazardous reactions.    |
| 10.4 Conditions to avoid:                | Not designated.                  |
| 10.5 Incompatible materials:             | Not designated.                  |
| 10.6 Hazardous decomposition products:   | Carbon monoxide, carbon dioxide. |

**11. Toxicological information:**

11.1 Information on toxicological effects

|   |   |
|---|---|
| When used appropriately no dangers to health are currently known              |   |
| Acute toxicity:   | No information available on systemic effects.   |
| Skin irritation:  | Danger of burns when handling hot and liquid waxes.   |
| Eye irritation:   | May harm when working with hot and liquid waxes.  |
| Sensitisation to the respiratory tract / skin:                                | Vapours have a very weak irritant effect.   |
| Germ cell mutagenicity:   | No information available  |
| Carcinogenicity:  | An extensive programme already realized in 1958 for clarifying the absence of carcinogenic properties in pure waxes was followed by the admission as food ingredient. |
| Reproductive toxicity:  | No information available  |
| Specific target organ toxicity (single exposure):                             | none  |
| Specific target organ toxicity (repeated exposure):                           | none  |
| Aspiration hazard:  | No information available  |
| Symptoms related to the physical, chemical and toxicological characteristics: |   |
| In case of ingestion:   | may have a laxative effect  |
| In case of skin contact:  | may cause burns when wax hot and liquid   |



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 5 of 6

Printing date: 19.05.2016

**Wax BioStar**

|                         |  |
|-------------------------|--|
| In case of inhalation:  | seems impossible   |
| In case of eye contact: | seek for medical advice when eye contact with hot and liquid wax |

**12. Ecological information:**

|   |   |
|---|---|
| 12.1 Toxicity:                            | No information available  |
| 12.2 Persistence and degradability:       | No data available. As to consistency and non-solubility a bioavailability is very unlikely. |
| 12.3 Bioaccumulative potential:           | Not significantly accumulative in living organisms.   |
| 12.4 Mobility in soil:                    | No data available.  |
| 12.5 Results of PBT- and vPvB assessment: | Not applicable.   |
| 12.6 Other adverse effects:               | None  |

**13. Disposal considerations:**

|  |   |
|--|---|
| 13.1 Waste treatment methods:                                | Contaminated packaging disposal<br>Recycle according to official regulations. |
| Waste codes / waste designations according to EWC / AVV:     | No waste codes applicable.  |
| Special precautions:   | None  |
| Other relevant provisions for waste treatment national or EC |   |

**14. Transport information:**

|   |   |
|---|---|
| 14.1 UN-no.   | Not applicable  |
| 14.2 UN proper shipping name<br>ADR/RID:<br>IMDG-Code / ICAO-TI / IATA-DGR:   | The product is not defined under ADR/RID regulations.<br>The product is not defined under IMDG regulations. |
| 14.3 Transport hazard class(es) / hazard labels:  | Not relevant.   |
| 14.4 Packing group:   | Not relevant.   |
| 14.5 Environmental hazards<br>Environmentally hazardous substance<br>marking<br>ADR/RID / IMDG-Code / ICAO-TI / IATA-DGR:                     | <input type="checkbox"/> yes / XXX no   |
| Marine Pollutant:   | <input type="checkbox"/> yes / XXX no   |
| 14.6 Special precautions for user:  | No additional information available.  |
| 14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code<br>Pollution category (X, Y or Z):<br>Type of ship (1, 2 or 3): | Not applicable.<br>Not applicable.  |

**15. Regulatory information:**

|  |  |
|--|--|
| 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture<br>EU-regulations, e.g.: |  |
| Directive /EC) no. 2037/2000 (substances that deplete the ozone layer):  | not listed                                 |
| Directive (EC) no. 850/2004 (persistent organic pollutants):   | not listed                                 |
| Directive (EC) no. 689/2008 (Import and export of dangerous chemicals):  | not listed                                 |
| Directive (EC) no. 648/2004 (Detergents Regulation):   | not listed                                 |
| National regulations, e.g.:  |  |
| Water hazard class (VwVwS) :   | not water pollutant / registration no. 268 |
| Technical instructions on Air Quality Control (TA-Luft):   | not applicable                             |



According to 1907/2006 (REACH)

Date of issue: 23.04.2015

Date of revision 22.01.16

Page 6 of 6

Printing date: 19.05.2016

**Wax BioStar**

|                                  |  |
|----------------------------------|--|
| Other relevant regulations:      | TRGS 500, 510 (Germany): storage class: 11 (flammable solids)      |
| 15.2 Chemical safety assessment: | For this product a chemical safety assessment was not carried out. |

**16. Other information:**

---

|  |   |
|--|---|
| Indication of changes:   | Regulation (EC) no. 1272/2008 (CLP, EU_GHS)<br>Section1: Emergency number, editorial revision   |
| Abbreviations and acronyms:  | n.a. Not applicable   |
| Key literature references and sources for data:  | Gestis Stoffdatenbank - <a href="http://gestis.itrust.de/">http://gestis.itrust.de/</a><br>Authority for chemicals <a href="http://www.baua.de">www.baua.de</a> |
| Classification for mixtures and used evaluation method according to regulation (EC) 1207/2008 [CLP]: | No classification   |
| Relevant R-, H- and EUH-phrases (number and full text):  | None of the above stated phrases is relevant.<br>P-phrases stated in section 2.2  |
| Training advice:   | Not applicable for qualified dental technicians.  |
| Further information:   | none  |

This information is based on our present state of knowledge. It should therefore not be construed as guaranteeing specific properties of the products described or their suitability for a particular application. We do not take any responsibility for possible sources of errors or sufficiency of this document. The user of the product accepts the sole responsibility to adhere to all herein mentioned regulations. It is the user's obligation to read this document carefully.